



ORTEC Software
File Structure Manual
for DOS and Windows® Systems

Advanced Measurement Technology, Inc.

a/k/a/ ORTEC®, a subsidiary of AMETEK®, Inc.

WARRANTY

ORTEC® DISCLAIMS ALL WARRANTIES OF ANY KIND, EITHER EXPRESSED OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE, NOT EXPRESSLY SET FORTH HEREIN. IN NO EVENT WILL ORTEC BE LIABLE FOR INDIRECT, INCIDENTAL, SPECIAL, OR CONSEQUENTIAL DAMAGES, INCLUDING LOST PROFITS OR LOST SAVINGS, EVEN IF ORTEC HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH DAMAGES RESULTING FROM THE USE OF THESE DATA.

Copyright © 2002, Advanced Measurement Technology, Inc. All rights reserved.

*ORTEC® is a registered trademark of Advanced Measurement Technology, Inc. All other trademarks used herein are the property of their respective owners.

TABLE OF CONTENTS

| | |
|---|----|
| 1. INTRODUCTION | 1 |
| 2. MESSAGE FILES | 3 |
| 2.1. DOS | 3 |
| 2.2. Windows | 3 |
| 3. EMULATOR FILES | 5 |
| 3.1. Integer Data Files (.CHN) | 5 |
| 3.1.1. Early Versions | 5 |
| 3.1.2. New Versions | 6 |
| 3.1.3. Old AlphaVision Version | 7 |
| 3.2. ROI Files | 8 |
| 3.3. Start Files | 8 |
| 3.4. Other Files | 8 |
| 4. ANALYSIS FILES | 9 |
| 4.1. Inform Structure | 9 |
| 4.2. Analysis Parameters Files | 9 |
| 4.2.1. Germanium Parameters | 9 |
| 4.2.2. Sodium Iodide Parameters | 13 |
| Analysis Parameters Extension 1 Record | 14 |
| Record 2 | 15 |
| Record 3 | 15 |
| 4.2.3. AlphaVision Parameters | 16 |
| 4.3. Calibration Files | 17 |
| 4.3.1. Calibration Description Record | 17 |
| 4.3.2. Calibration Parameter Records | 17 |
| 4.3.2.1. Germanium Calibration | 17 |
| First Certificate Record | 19 |
| First Total Record | 20 |
| Second Calibration Data Record | 20 |
| 4.3.2.2. Sodium Iodide Calibration | 22 |
| 4.4. Detector Description Files | 23 |
| 4.4.1. Detector Description Record | 23 |
| 4.5. Old Geometry Correction Files | 23 |
| 4.5.1. Geometry Correction Data Record | 24 |
| 4.5.2. Geometry Correction Description Record | 24 |
| 4.6. New Geometry Correction Files | 25 |
| 4.7. Sample Description Files | 25 |
| 4.7.1. Sample Description Record | 26 |

| | |
|--|----|
| 4.8. ROI File | 26 |
| 4.8.1. Record 1 | 26 |
| 4.8.2. ROI Data Record 1 and 2 | 26 |
| 4.9. Old Absorption Correction Files | 27 |
| 4.9.1. Absorption Correction Data Record | 28 |
| 4.9.2. Absorption Correction Description Record | 28 |
| 4.10. New Absorption Correction Files | 29 |
| 4.11. Real Format Spectrum Files | 29 |
| 4.11.1. Record 1 | 29 |
| 4.11.2. Acquisition Information Record | 31 |
| 4.11.3. Sample Description Record | 31 |
| 4.11.4. Detector Description Record | 32 |
| 4.11.5. First, Second, and Third Analysis Parameter | 32 |
| 4.11.6. Absorption Correction Description Record | 32 |
| 4.11.7. Absorption Correction Data Record | 32 |
| 4.11.8. Geometry Correction Description Record | 33 |
| 4.11.9. Geometry Correction Data Record | 33 |
| 4.11.10. Calibration Description Record, Calibration Data Records, and Efficiency Pairs Records | 33 |
| 4.11.11. Region-of-Interest Records | 33 |
| Record 1 | 34 |
| Record 2 | 34 |
| Record 3 | 34 |
| 4.11.12. Hardware Parameters Records | 35 |
| Record 1 | 35 |
| Record 2 | 36 |
| Record 3 | 37 |
| Record 4 | 38 |
| Multi-Nuclide MDA Preset Storage | 38 |
| 4.11.13. Personality Record | 39 |
| 4.11.14. Table Description Record | 40 |
| 4.11.15. Table Data Record | 40 |
| 4.11.16. Spectrum Records | 40 |
| 4.12. Integer Format Spectrum Files | 41 |
| Record 1 | 41 |
| 4.12.1. Other Records | 42 |
| 4.12.2. Spectrum Records | 42 |
| 4.13. Net Spectrum Files | 42 |
| Record 1 | 43 |
| 4.14. AlphaVision Spectrum Files | 44 |
| 4.15. Germanium Unformatted Output Files | 45 |

| | |
|--|----|
| 4.15.1. Records 1 and 2 | 46 |
| 4.15.1.1. Old UFO Records 1 and 2 | 46 |
| Record 1 | 46 |
| Record 2 | 47 |
| Other Records | 48 |
| 4.15.1.2. New UFO Records 1 and 2 | 48 |
| Record 1 | 48 |
| Record 2 | 50 |
| Other Records | 51 |
| 4.15.2. UFO File CSI Records | 51 |
| Records 1 and 3 | 51 |
| Records 2 and 4 | 52 |
| 4.15.3. UFO File GEN Record | 52 |
| 4.15.4. UFO Absorption Description Record | 53 |
| 4.15.5. UFO Geometry Description Record | 54 |
| 4.15.6. UFO File Nuclide Records | 55 |
| 4.15.6.1. Old UFO File Nuclide Records | 55 |
| 4.15.6.2. New UFO File Nuclide Records | 56 |
| 4.15.7. UFO File Peak Records | 58 |
| 4.15.7.1. Old UFO File Peak Records | 58 |
| 4.15.7.2. New UFO File Peak Records | 60 |
| 4.16. Sodium Iodide Unformatted Output (UFO) File Peak Records | 63 |
| 4.16.1. UFO Records 1 and 2 | 63 |
| Record 1 | 63 |
| Record 2 | 64 |
| 4.16.2. UFO File CSI Records | 65 |
| Records 1 and 3 | 65 |
| Records 2 and 4 | 66 |
| 4.16.3. UFO File GEN Record | 66 |
| 4.16.4. UFO File Nuclide Records | 67 |
| 4.16.5. UFO File Peak Records | 68 |
| 4.17. AlphaVision Output File | 70 |
| 5. ASCII FORMAT FILES | 73 |
| 5.1. SPE Spectrum Files | 73 |
| 5.2. IEEE Spectrum Format | 73 |
| 6. LIBRARY FILES | 75 |
| 6.1. Library Files (Germanium) | 75 |
| 6.1.1. Header Record | 75 |
| 6.1.2. Nuclide Records | 76 |

| | |
|--|----|
| 6.1.3. Peak Records | 77 |
| 6.2. PBC Files (Germanium) | 77 |
| 6.2.1. Header Record | 78 |
| 6.2.2. Nuclide Records | 79 |
| 6.2.3. Peak Records | 79 |
| 6.3. Library Files in NuclideNavigator | 80 |
| 6.4. Library Files in AlphaVision | 80 |
| 6.4.1. Nuclide Records | 80 |
| 7. DOS SYSTEM FILES | 83 |
| 7.1. Sequence Number File | 83 |
| 7.2. Configuration File | 83 |
| 7.2.1. Configuration Records | 84 |
| 7.3. Report File | 84 |
| 8. TABLE FILES | 85 |
| 8.1. Header Record | 85 |
| 8.2. Title Records | 85 |
| 8.3. Nuclide Records | 85 |
| 8.3.1. PBC Type | 86 |

1. INTRODUCTION

This manual describes the file structure for the files used in the ORTEC MCA control software and the analysis software. These files are used to transfer information among the programs in these products and to store data for future use. Developers who wish to design programs to use these files or to create files that will be used by ORTEC software must adhere to the specifications given.

There are six types of files used:

- Message files (**.MSG**)
- Integer spectrum data files (**.CHN**)
- Analysis files (various extensions)
- Library files (**.LIB**)
- System files (various extensions)
- Table files (**.TBL**)

Other file types (**.BAT**, **.EXE**, **.COM**) are system files and defined in other manuals.

All real numbers used in these files conform to the IEEE real data type as implemented by Microsoft® and Intel®.

The files described here, especially the spectrum files and the analysis results files, have been expanded to include more information about the details of the spectrum collection and the analysis options. This was done to enable customers to better maintain their file archives and to support their QA requirements.

Programs written for the DOS operating system are no longer supported because they do not pass the Y2K tests. References to DOS in this manual are for the convenience of those using old programs and does not imply support.

2. MESSAGE FILES

2.1. DOS

The message files contain the operator prompts, keyword replies, and error messages for the programs. There is no internal file header. They are created as 64-byte fixed-length direct access records (i.e., there are no interrecord delimiters). Program MSGCOM is used to compile the random length sequential access records (i.e., there are interrecord delimiters), as from a text editor, into this format. The extension is **.MSG**.

2.2. Windows

The message (**.TXT**) and initialization files (**.INI**) contain the operator prompts, keyword replies, startup conditions, and error messages for the programs. There is no internal file header. They are created as random length, sequential access records (i.e., there are interrecord delimiters). They are read by the program as text files and require no processing.

3. EMULATOR FILES

3.1. Integer Data Files (.CHN)

.CHN files can be produced by the pulse-height analysis (PHA) programs (A63-BI, A64-BI, A65-BI, and A65-B32), Transfer Program (A48-BI), Applications Manager (A18-BI), analysis programs (A30-BI, A34-BI, A66-BI, and A66-B32) and others.

The .CHN integer data files contain the channel-by-channel contents of the MCB. The header is 32-bytes long and contains the following:

| <u>Byte Offset</u> | <u>Byte Length</u> | <u>Use</u> |
|------------------------|------------------------|--|
| 0 | 2 | Must be -1 |
| 2 | 2 | MCA number or Detector number |
| 4 | 2 | Segment number (set to 1 in UMCBI) |
| 6 | 2 | ASCII seconds of start time |
| 8 | 4 | Real Time (increments of 20 ms) (4-byte integer) |
| 12 | 4 | Live Time (increments of 20 ms) (4-byte integer) |
| 16 | 8 | Start date as ASCII DDMMYY* or binary zeros, if not known. The * character should be ignored if it is not a "1". If it is a "1", it indicates the data is after the year 2000. |
| 24 | 4 | Start time as ASCII HHMM or binary zeros, if not known (see Byte 6 above) |
| 28 | 2 | Channel offset of data |
| 30 | 2 | Number of channels (length of data) |

The next part of the file contains the spectrum stored as 4-byte integers. There are no record separators in the file. The number of spectrum records is determined by the number of channels in the spectrum.

3.1.1. Early Versions

The last part of the spectrum file contains additional descriptive information about the spectrum, as follows (the byte offsets are relative to the end of the spectrum) for all DOS versions, for GammaVision® V2 and below, MAESTRO® for Windows V3 and below:

| <u>Byte Offset</u> | <u>Byte Length</u> | <u>Use</u> |
|------------------------|------------------------|--|
| 0 | 2 | Must be -101 |
| 2 | 2 | Reserved |
| 4 | 4 | Energy calibration zero intercept, 0.0 for uncalibrated spectrum (4-byte real) |

| | | |
|-----|-----|--|
| 8 | 4 | Energy calibration slope, 1.0 for uncalibrated spectrum (4-byte real) |
| 12 | 4 | Reserved |
| 16 | 4 | Peak shape calibration zero intercept, 1.0 for uncalibrated spectrum (4-byte real) |
| 20 | 4 | Peak shape calibration slope, 0.0 for uncalibrated spectrum (4-byte real) |
| 24 | 232 | Reserved |
| 256 | 1 | Length of detector description |
| 257 | 63 | Detector description |
| 320 | 1 | Length of sample description |
| 321 | 63 | Sample description |
| 384 | 128 | Reserved (The total length is 512 bytes) |

This file is also created by STORE (A18-BI) and XFER (A48-BI).

3.1.2. New Versions

The last part of the spectrum file contains additional descriptive information about the spectrum, as follows (the byte offsets are relative to the end of the spectrum) for GammaVision V2.2 and later, MAESTRO V3.1 and later, MicroMCB and all versions of programs for Windows 95/98/NT:

| <u>Byte Offset</u> | <u>Byte Length</u> | <u>Use</u> |
|------------------------|------------------------|--|
| 0 | 2 | Must be -102 |
| 2 | 2 | Reserved |
| 4 | 4 | Energy calibration zero intercept, 0.0 for uncalibrated spectrum (4-byte real) |
| 8 | 4 | Energy calibration slope, 1.0 for uncalibrated spectrum (4-byte real) |
| 12 | 4 | Energy calibration quadratic term, 0.0 for uncalibrated spectrum (4-byte real) |
| 16 | 4 | Peak shape calibration zero intercept, 1.0 for uncalibrated spectrum (4-byte real) |
| 20 | 4 | Peak shape calibration slope, 0.0 for uncalibrated spectrum (4-byte real) |
| 24 | 4 | Peak shape calibration quadratic term, 0.0 for uncalibrated spectrum (4-byte real) |
| 28 | 228 | Reserved |
| 256 | 1 | Length of detector description |
| 257 | 63 | Detector description |
| 320 | 1 | Length of sample description |
| 321 | 63 | Sample description |
| 384 | 128 | Reserved (The total length is 512 bytes) |

3.1.3. Old AlphaVision Version

The old (A36-BI) versions of AlphaVision used the following **.CHN** files. This file type is not used in the new (A36-B32) versions. See Section 4.14 for the new AlphaVision spectrum files.

| Byte Offset | Byte Length | Use |
|----------------|----------------|--|
| 0 | 2 | Must be -102 |
| 2 | 2 | Reserved |
| 4 | 4 | Energy calibration zero intercept, 0.0 for uncalibrated spectrum (4-byte real) |
| 8 | 4 | Energy calibration slope, 1.0 for uncalibrated spectrum (4-byte real) |
| 12 | 4 | Energy calibration quadratic term, 0.0 for uncalibrated spectrum (4-byte real) |
| 16 | 4 | Peak shape calibration zero intercept, 1.0 for uncalibrated spectrum (4-byte real) |
| 20 | 4 | Peak shape calibration slope, 0.0 for uncalibrated spectrum (4-byte real) |
| 24 | 4 | Peak shape calibration quadratic term, 0.0 for uncalibrated spectrum (4-byte real) |
| 28 | 228 | Reserved |
| 256 | 1 | Length of detector description |
| 257 | 63 | Detector description |
| 320 | 1 | Length of sample description |
| 321 | 63 | Sample description |
| 384 | 4 | AlphaVision ID must be 0x53495641 |
| 388 | 32 | Sample type name |
| 420 | 10 | ASCII collection date |
| 430 | 10 | ASCII collection time |
| 440 | 4 | Total volume of sample |
| 444 | 4 | Aliquot volume of sample |
| 448 | 4 | Tracer amount (DPM) |
| 450 | 2 | Volume units |
| 454 | 4 | Detector efficiency |
| 458 | 4 | Old calibration intercept |
| 462 | 4 | Old calibration slope |
| 466 | 4 | Old calibration shape |
| 470 | 4 | Old calibration efficiency |
| 474 | 4 | Old background counts |
| 478 | 4 | Old background CPM |
| 482 | 30 | Group Name |

3.2. ROI Files

This file is created by SAVE ROI and contains a list of the start and stop channels for the ROIs in the display. The file contents are as follows:

| <u>Byte Offset</u> | <u>Byte Length</u> | <u>Use</u> |
|------------------------|------------------------|--------------------------------------|
| 0 | 2 | Must be -2 |
| 2 | 2 | Start channel number of first ROI |
| 4 | 2 | Stop channel number of first ROI |
| . | . | . |
| . | . | . |
| . | . | . |
| | | Continue for all ROIs in the display |
| n | 2 | Start channel = -1 is end of data |

3.3. Start Files

This file is written by the DOS Emulator (A63-BI) to store the start time and date of the MCB/Segment data. It is used by Emulator when storing the data on disk. There is no header block. This file is not used in any Windows program. The contents are:

| <u>Byte Offset</u> | <u>Type</u> | <u>Use</u> |
|------------------------|-------------|-------------------------|
| 0 | I*2 | File type = -4 |
| 2 | | Unused |
| 16 | C*8 | Start date as MMM-DD-YY |
| 24 | C*8 | Start time as HH:MM:SS |
| 28 | | Unused |

This file is put in the default directory.

3.4. Other Files

The file **PRECAL.MCB** is used for internal data transfer in the DOS Emulator. The contents are reserved.

4. ANALYSIS FILES

4.1. Inform Structure

The analysis files are written in a defined format called the inform structure. This is a format that allows great flexibility in the use and contents of the file. The header record contains the definition of the contents and their position(s) in the file. As all of the inform files are dynamically allocated, the first record must be carefully used, and must be modified if the file is modified. If the particular record type doesn't exist in the file, the pointer is zero or negative. The recommended (or default) extension is given in the following. The file extension does not, however, override the file-type bits set in the file. All inform files have a record length of 128 bytes with no interrecord delimiter.

Significant changes have been made for analysis software files for Windows 95/98/NT.

4.2. Analysis Parameters Files

An analysis parameters file is created by GERPAR (B30) or MINPAR (A34) and contains the analysis parameters used in the analysis. The file extension is **.PRM**. The first record in this file contains the pointers to the four analysis parameters records (Words 9 through 12).

| Word Number | Use |
|----------------|---|
| 1 | Must be 1 |
| 2 | Must be 32 |
| 3–8 | Reserved |
| 9 | First analysis parameters record pointer |
| 10 | Second analysis parameters record pointer |
| 11 | Third analysis parameters record pointer |
| 12 | Hardware parameters record pointer |
| 13–27 | Reserved |
| 28 | Maximum record number ever used |
| 29 | Maximum record number in use |
| 30–69 | Reserved |

4.2.1. Germanium Parameters

The first two records are used for germanium analysis (A34 and B30). The files are not produced by GammaVision but the analysis parameters records in the spectrum and other files are in this format. The first germanium record contains the following data:

| <u>Word</u> <u>Number</u> | <u>Local</u> <u>Name</u> | <u>Type</u> | <u>Use</u> | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
|------------------------------|-----------------------------------|-------------|--|------------|--------------------|-------|-------------------------------|-------|-----------------------------------|-------|---------------------------|-------|--------------------|-------|------------|-------|----------------------|-------|-------------------|-------|----------------------|-------|-----------------------|-------|----------------------|--------|---------------------|--------|--------------------------|-----------|----------|
| 1 | GELI | I*2 | Bit 0 must be 1 for germanium. Bit 14 = 1 for extended analysis records. | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 2 | FCHAN | I*2 | Start channel for analysis | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 3 | LCHAN | I*2 | Stop channel for analysis | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 4 | SIGMA | I*2 | Standard deviation for uncertainty (1, 2, or 3) | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 5 | UNPEKS | I*2 | Maximum number of unused peaks allowed | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 6 | ENGCAL | I*2 | True if spectrum is energy-calibrated | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 7 | DECDUR | L*2 | Decay during acquisition flag | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 8 | FORM | I*2 | 1 = Count as %, 2 = Count as activity, 3 = Count and total as %, 4 = Count and total as activity | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 9 | GAMMA3 | L*2 | (Old) True to enable peak interference correction (Gamma 3) | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | I*2 | (New) Analysis flag controls as follows (true if bit set to 1) | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | <table><tr><th><u>Bit</u></th><th><u>Description</u></th></tr><tr><td>Bit 0</td><td>Enable library peak stripping</td></tr><tr><td>Bit 1</td><td>Enable average energy calculation</td></tr><tr><td>Bit 2</td><td>Enable iodine equivalence</td></tr><tr><td>Bit 3</td><td>Write ASCII report</td></tr><tr><td>Bit 4</td><td>Report RPG</td></tr><tr><td>Bit 5</td><td>Run at high priority</td></tr><tr><td>Bit 6</td><td>Directed fit flag</td></tr><tr><td>Bit 7</td><td>Use GV Report Writer</td></tr><tr><td>Bit 8</td><td>Manual peak stripping</td></tr><tr><td>Bit 9</td><td>Use collection dates</td></tr><tr><td>Bit 10</td><td>Enable TCC analysis</td></tr><tr><td>Bit 11</td><td>Display analysis results</td></tr><tr><td>Bit 12–16</td><td>Reserved</td></tr></table> | <u>Bit</u> | <u>Description</u> | Bit 0 | Enable library peak stripping | Bit 1 | Enable average energy calculation | Bit 2 | Enable iodine equivalence | Bit 3 | Write ASCII report | Bit 4 | Report RPG | Bit 5 | Run at high priority | Bit 6 | Directed fit flag | Bit 7 | Use GV Report Writer | Bit 8 | Manual peak stripping | Bit 9 | Use collection dates | Bit 10 | Enable TCC analysis | Bit 11 | Display analysis results | Bit 12–16 | Reserved |
| <u>Bit</u> | <u>Description</u> | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Bit 0 | Enable library peak stripping | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Bit 1 | Enable average energy calculation | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Bit 2 | Enable iodine equivalence | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Bit 3 | Write ASCII report | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Bit 4 | Report RPG | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Bit 5 | Run at high priority | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Bit 6 | Directed fit flag | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Bit 7 | Use GV Report Writer | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Bit 8 | Manual peak stripping | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Bit 9 | Use collection dates | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Bit 10 | Enable TCC analysis | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Bit 11 | Display analysis results | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Bit 12–16 | Reserved | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 10 | MPCON | L*2 | Maximum permissible concentration on/off | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 11 | PBCOR | L*2 | Peaked background correction on/off | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 12 | RNDSUM | L*2 | Random summing on/off | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 13 | DECAY | L*2 | Decay correction on/off | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 14 | GEOM | L*2 | Geometry correction on/off | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 15 | FACCOR | R*4 | Reserved | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 17 | EFCALB | L*2 | True if spectrum is efficiency-calibrated | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 18 | MDATYP | I*2 | MDA type, 1 to 20 (see B30 manual) | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 19 | LIBNOT | L*2 | True if analysis library found | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 20 | LSTISO | L*2 | True if isotopic abundance output requested | | | | | | | | | | | | | | | | | | | | | | | | | | | | |

| | | | |
|----|--------|-----|---|
| 21 | FIVE | I*2 | 1 = Automatic background, 2 = five-point background width only, 3 = 3-point background, 4 = minimum point |
| 22 | SENS | R*4 | Sensitivity |
| 24 | RECFCT | R*4 | Reserved |
| 25 | DETIDN | I*2 | Detector identification number |
| 26 | MATISO | L*2 | True if isotopic/peak output requested |
| 27 | MATCHN | I*2 | Output print control, true if bit set to 1. |
| | | | <u>Bit</u> <u>Description</u> |
| | | | Bit 0 Print library peak list |
| | | | Bit 1 Print unknown peak list |
| | | | Bit 2–16 Reserved |
| 28 | RSFCT | R*4 | Random summing factor |
| 30 | SORFCT | R*4 | Absorption factor |
| 32 | DAYS | R*8 | Decay correction date in DECDAY format. Days since 1Jan79. The example time date, 2Jan79 12:00:00, is 1.5 |
| 36 | SAMWGH | R*4 | Sample weight |
| 38 | UFCTN | R*4 | Units factor numerator |
| 40 | UFCTD | R*4 | Units factor denominator |
| 42 | IRRAD | R*4 | Reserved |
| 44 | PKRPRO | R*4 | Additional random error |
| 46 | GERPRO | R*4 | Additional systematic error |
| 48 | UNITS | C*1 | Units name in ASCII (14 characters) |
| 55 | CORDAT | C*1 | Date in ASCII (DD-MMM-YY*) (10 characters) |
| 60 | CORTIM | C*1 | Time in ASCII (HH:MM:SS) (8 characters) |
| 64 | SORP | L*2 | Absorption correction flag |

The DECDAY format is defined to be a double precision real number that represents the number of days since 1-JAN-79. For dates before this date the number is negative.

The extended analysis record follows the above and contains:

| <u>Word</u> | <u>Use</u> |
|---------------|---------------------------------|
| <u>Number</u> | |
| 1–2 | Fraction Limit test in percent |
| 3–4 | Unknown peak cutoff (REN-B32-G) |
| 5 | Divide activity by weight |
| 6–7 | Half-lives decay limit |
| 8–9 | Activity range factor |

| | |
|-------|--------------------------------------|
| 10–11 | Minimum step background energy |
| 12 | Second MDA type |
| 13–16 | Analysis version |
| 17–32 | Reserved |
| 33–48 | Second library for manual peak strip |
| 49–64 | Third library for manual peak strip |

The second germanium record contains the library names used in the analysis.

| <u>Word</u> <u>Number</u> | <u>Local</u> <u>Name</u> | <u>Type</u> | <u>Use</u> |
|------------------------------|-----------------------------|-------------|---------------------------------------|
| 1 | LIB1 | C*32 | First nuclide library filename |
| 17 | LIB2 | C*32 | Attenuation filename |
| 33 | LIB3 | C*32 | Geometry correction filename |
| 49 | PBCFIL | C*32 | Peaked background correction filename |

If long filenames are enabled, then the following 8 sequential records are defined as follows:

| | |
|-----|---|
| 1+2 | 256 characters for LIB1 |
| 3+4 | 256 characters for Attenuation filename |
| 5+6 | 256 characters for Geometry filename |
| 7+8 | 256 characters for PBC filename |

The third germanium record contains more filenames used in the analysis.

| <u>Word</u> <u>Number</u> | <u>Local</u> <u>Name</u> | <u>Type</u> | <u>Use</u> |
|------------------------------|-----------------------------|-------------|--------------------------|
| 1 | RPGPROG | C*32 | MPCFIL |
| 17 | RPTOUT | C*32 | Report filename |
| 33 | SUSLIB | C*32 | Suspected nuclide file |
| 49 | ANLPROG | C*32 | Analysis executable name |

If long filenames are enabled, then the following 8 sequential records are defined as follows:

| | |
|-----|--------------------------|
| 1+2 | MPCFIL filename |
| 3+4 | Report filename |
| 5+6 | Suspected nuclide file |
| 7+8 | Analysis executable name |

4.2.2. Sodium Iodide Parameters

Analysis parameters record structure:

| <u>Word</u> <u>Number</u> | <u>Local</u> <u>Name</u> | <u>Type</u> |
|------------------------------|-----------------------------|--|
| 1 | GELI | Low byte must be 2; if bit 14 is 1, the extended analysis record is included |
| 2 | FCHAN | Start channel for analysis |
| 3 | LCHAN | Last channel |
| 4 | SIGMA | Standard deviation for uncertainty (1, 2, or 3) |
| 5 | UNPEKS | Maximum number of unused peaks allowed |
| 6 | ENGCAL | Spectrum is energy calibrated? |
| 7 | DECDUR | Decay during acquisition? |
| 8 | FORM | (1 = Count as %, 2 = Count as activity, 3 = Count and total %, 4 = Count and total activity) |
| 9 | GAMMA3 | |
| | | <u>Bit</u> <u>Description if = 1</u> |
| | | 0 Enable library peak stripping |
| | | 1 Enable average energy calc |
| | | 2 Enable iodine equivalence |
| | | 3 Save report as file |
| | | 4 Report RPG |
| | | 5 Run at high priority |
| | | 6 Directed fit flag |
| | | 7 Use the GammaVision Report Writer |
| | | 8 Fix the linear background fit |
| | | 9 Report to printer |
| | | 10 Report MDA (RENAISSANCE) |
| | | 11 Report individual detectors (RENAISSANCE-32) |
| 10 | MPCON | Maximum permissible concentration |
| 11 | PBCOR | Peaked background correction |
| 12 | RNDSUM | Random summing correction |
| 13–14 | DECAY | Calculation of ICRFAC |
| 15–16 | FACCOR | (Usually 0.5) acceptance width for identification of peak to library |
| 17 | EFFCALB | Spectrum is efficiency calibrated |
| 18 | MDATYP | MDA type |
| 19 | LIBNOT | Analysis library found? |

| | | |
|-------|---------|---|
| 20 | LSTISO | Isotopic abundance listing requested |
| 21 | FIVE: | Background type |
| 22-23 | SENS | Peak error cut off |
| 24-25 | | Peak cutoff |
| 26 | MATISO | Isotopic/peak output requested? |
| 27 | MATCHN | Energy/peak output requested |
| 28-29 | | Linear background energy |
| 28-29 | RSFCT | Random summing factor |
| 30-31 | SORFCT | Absorption factor |
| 32-35 | DAYS | DECDAY format of decay correction date |
| | | COR DAT,CORTIM |
| 36 | SAMWGH | Sample weight or volume |
| 38 | UFCTN | Units factor numerator |
| 40 | UFCTD | Units factor denominator |
| 42-43 | | Sample collection time |
| 44-45 | | Fraction limit |
| 46-47 | | Decay limit |
| 48-54 | UNITS | Units name |
| 55-59 | COR DAT | Decay correction date (DD-MMM-YY*) (sample collection date) |
| 60-63 | CORTIM | Decay correction time (HH:MM:SS) (sample collection time) |
| 64 | SORP | Absorption correction |

Analysis Parameters Extension 1 Record

This record follows next if bit 14, word 1 above is set to 1.

| Word Number | Use |
|----------------|--|
| 1 | 1 = use peak cutoff, 0 = critical level cutoff |
| 2 | Directed fit flag |
| 3-4 | Fraction limit test in percent |
| 5 | Divide activity by weight (REN-B32 only) |
| 5-64 | Reserved |

Record 2

The second sodium iodide record contains the library names used in the analysis.

| <u>Word Number</u> | <u>Local Name</u> | <u>Type</u> | <u>Use</u> |
|------------------------|-----------------------|-------------|---------------------------------------|
| 1 | LIB1 | C*32 | First nuclide library filename |
| 17 | LIB2 | C*32 | Attenuation filename |
| 33 | LIB3 | C*32 | Geometry correction filename |
| 49 | PBCFIL | C*32 | Peaked background correction filename |

If long filenames are enabled, then the following 8 sequential records are defined as follows:

| | |
|-----|---|
| 1+2 | 256 characters for LIB1 |
| 3+4 | 256 characters for attenuation filename |
| 5+6 | 256 characters for geometry filename |
| 7+8 | 256 characters for PBC filename |

Record 3

The third sodium iodide record contains more filenames used in the analysis.

| <u>Word Number</u> | <u>Local Name</u> | <u>Type</u> | <u>Use</u> |
|------------------------|-----------------------|-------------|--------------------------|
| 1 | RPGPROG | C*32 | MPCFIL |
| 17 | RPTOUT | C*32 | Report filename |
| 33 | SUSLIB | C*32 | Suspected nuclide file |
| 49 | ANLPROG | C*32 | Analysis executable name |

If long filenames are enabled, then the following 8 sequential records are defined as follows:

| | |
|-----|--------------------------|
| 1+2 | MPCFIL filename |
| 3+4 | Report filename |
| 5+6 | Suspected nuclide file |
| 7+8 | Analysis executable name |

4.2.3. AlphaVision Parameters

AlphaVision analysis record 1.

| <u>Word</u> <u>Number</u> | <u>Use</u> |
|------------------------------|------------------------------------|
| 0 | must be 4 |
| 1 | AlphaVision ID, must be Ox53495641 |
| 2 | Sample type name |
| 18 | ASCII collection date |
| 23 | ASCII collection time |
| 27 | Total Volume of sample |
| 29 | Aliquot Volume of sample |
| 31 | Tracer amount (DPM) |
| 33 | Tracer units |
| 34 | Detector Efficiency |
| 36 | Old calibration intercep |
| 38 | Old calibration slope |
| 40 | Old calibration shape |
| 42 | Old calibration efficiency |
| 44 | Old background counts |
| 46 | Old background CPM |
| 48 | Group name |
| 54–255 | Reserved |
| 256 | Calibration intercept |
| 258 | Calibration slope (1.0 = none) |
| 260 | Calibration shape (5.93) |
| 262 | Peak shape intercept (1.0 = none) |
| 264 | Peak shape slope (0.0 = none) |
| 266 | Peak shape |
| 268 | Calibration units |
| 270 | Efficiency |
| 272–329 | Reserved |
| 330 | Detector description length |
| 331–362 | Detector description |
| 363 | Sample description length |
| 364–385 | Sample description |
| 386–511 | Reserved |

Analysis records 2 and 3 are not used in AlphaVision.

4.3. Calibration Files

A calibration file is created by the calibration program and contains the calibration information used in the analysis programs. The filename extension should be **.CLB**. This file contains the calibration description (if any), the calibration parameter records, and the efficiency pairs records (if any). The calibration file contents are:

| <u>Word</u> <u>Number</u> | <u>Use</u> |
|------------------------------|---|
| 1 | Must be 1 |
| 2 | Must be 256 |
| 3 | Contents flag Bit 0 = 1 for long filename Bit 1 = 1 for ZDT spectrum and ROI in file. |
| 16 | Reserved |
| 17 | Calibration description record pointer |
| 18 | First calibration data record pointer |
| 19 | Second calibration data record pointer |
| 20 | Efficiency pairs record pointer (first record) |
| 21 | Reserved |
| 22 | Energy pairs record pointer (first record) |
| 23 | Number of energy pairs records (first record) |
| 24–27 | Reserved |
| 28 | Maximum record number ever used |
| 29 | Maximum record number in use |
| 30 | Number of efficiency pairs records (See Word 20) |
| 31–64 | Reserved |

4.3.1. Calibration Description Record

The calibration description record contains the user-supplied description of the calibration parameters. The record is divided into two 64-character lines of ASCII text. The first 64 characters are used for the energy calibration and the second 64 characters for the efficiency calibration.

4.3.2. Calibration Parameter Records

The first calibration data record contains the various coefficients produced by the calibration program, as described below. Note that the FWHM coefficients (FC(i)) produce the FWHM in channels, while the programs usually display the FWHM in energy units.

4.3.2.1. Germanium Calibration

| <u>Word Number</u> | <u>Local Name</u> | <u>Type</u> | <u>Use</u> |
|------------------------|-----------------------|-------------|---|
| 1 | AFIT | I*2 | Above knee efficiency calibration fit type |
| 2 | BFIT | I*2 | Below knee efficiency calibration fit type |
| 3 | EFFPRS | I*2 | Number of efficiency pairs |
| 4 | NCH | I*2 | Number of channels in spectrum |
| 5 | KNEE | R*4 | Detector knee (keV) |
| 7 | ASIG | R*4 | 2-sigma uncertainty above knee |
| 9 | BSIG | R*4 | 2-sigma uncertainty below knee |
| 11 | EC(1) | R*4 | Energy vs. channel coefficient A |
| 13 | EC(2) | R*4 | Energy vs. channel coefficient B |
| 15 | EC(3) | R*4 | Energy vs. channel coefficient C |
| 17 | FC(1) | R*4 | FWHM vs. channel coefficient A |
| 19 | FC(2) | R*4 | FWHM vs. channel coefficient B |
| 21 | FC(3) | R*4 | FWHM vs. channel coefficient C |
| 23 | PE(1) | R*4 | Above knee efficiency vs. energy coefficient A or polynomial coefficient (1) |
| 25 | PE(2) | R*4 | Above knee efficiency vs. energy coefficient B or polynomial coefficient (2) |
| 27 | PE(3) | R*4 | Above knee efficiency vs. energy coefficient C or polynomial coefficient (3) |
| 29 | SE(1) | R*4 | Below knee efficiency vs. energy coefficient A or polynomial coefficient (4) |
| 31 | SE(2) | R*4 | Below knee efficiency vs. energy coefficient B or polynomial coefficient (5) |
| 33 | SE(3) | R*4 | Below knee efficiency vs. energy coefficient C or polynomial coefficient (6) |
| 35 | FWHTYP | I*2 | FWHM type |
| 36 | PETYPE | | True for p-type |
| 37 | | | MAESTRO peak-search sensitivity |
| 38 | ENGPRS | I*2 | Number of energy pairs |
| 39 | DETNUM | I*2 | Detector number |
| 40 | NBKNEE | I*2 | Number of calibration points below knee |
| 41 | ENA2 | R*4 | Temp energy calibration |
| 43 | ENB2 | R*4 | Temp energy calibration |
| 45 | ENC2 | R*4 | Temp energy calibration |
| 47 | CALUNC | R*4 | Calibration source uncertainty |
| 49 | CALDIF | R*4 | Energy calibration difference |
| 51 | R(7) | R*4 | Polynomial coefficient 7 |
| 53 | R(8) | R*4 | Polynomial coefficient 8 |

| | | | |
|-------|-----------|-----|-----------------------------------|
| 55 | R(9) | R*4 | Polynomial coefficient 9 |
| 57 | R(10) | R*4 | Polynomial coefficient 10 |
| 59–60 | | | Low channel FWHM error |
| 61–62 | | | High channel FWHM error |
| 63 | | | Low channel limit for calibrating |
| 64 | STYPEFLAG | I*2 | True = next record has TCC data |

This record contains the TCC calibration information and the standard efficiency certificate data

| <u>Word</u> | <u>Use</u> |
|-------------|--|
| 1–2 | Number of TCC calibration records including this one |
| 3 | First certificate record pointer |
| 4 | First total record pointer |
| 5–6 | Reserved |
| 7–8 | Number of certificate energies |
| 9–10 | Number peak-to-total pairs |
| 11–12 | Reserved |
| 13–16 | Peak-to-total fit intercept |
| 17–20 | Peak-to-total fit slope |
| 21–24 | LTS efficiency intercept |
| 25–28 | LTS efficiency slope |
| 29–32 | LTS efficiency quadratic |
| 33–34 | Live-time preset for energy calibration |
| 35–36 | Live-time preset for efficiency calibration |
| 37–38 | Live-time preset for background |
| 39–40 | Calibration SPC live time |
| 41 | Absorber flag: 0 = without absorber 1 = with absorber |
| 42 | TCC calibration method 0 = multiple point source 1 = single point source 2 = single extended source |
| 43–47 | TCC calibration date |
| 48–51 | TCC calibration time |
| 52–64 | Unused |

First Certificate Record

| Word | |
|---------------|--------------------------------------|
| <u>Number</u> | <u>Use</u> |
| 1–4 | Isotope name |
| 5–8 | Isotope half-life in days |
| 9–12 | Peak energy |
| 13–16 | Branching ratio (%) |
| 17–20 | Measured value |
| 21–24 | Fitted value |
| 25–28 | Error in fit |
| 29–33 | Certificate date |
| 34–37 | Certificate time |
| 38 | Activity units |
| 39–40 | Live time for peak-to-total analysis |
| 41–64 | Unused |

First Total Record

The total record contains the peak-to-total ratio for the isotope and energy given. There are two entries per record. The record is repeated for as many times as needed for the complete list.

| Word | |
|---------------|--|
| <u>Number</u> | <u>Use</u> |
| | <i>Words 1–32 are for the first isotope in record</i> |
| 1– 4 | Isotope name (ASCII) |
| 5–8 | Peak energy |
| 9–12 | Fitted value |
| 13–16 | Peak-to-total ratio |
| 17–20 | Difference between the fit and ratio |
| 21–22 | Live time |
| 23–32 | Unused |
| | <i>Words 33–64 are for the second isotope in record</i> |
| 33–36 | Isotope name (ASCII) |
| 37–40 | Peak energy |
| 41–44 | Fitted value |
| 45–48 | Peak-to-total ratio |
| 49–52 | Difference between the fit and ratio |
| 53–54 | Live time |
| 55–64 | Unused |

Second Calibration Data Record

| <u>Word</u> <u>Number</u> | <u>Local</u> <u>Name</u> | <u>Type</u> | <u>Use</u> |
|------------------------------|-----------------------------|-------------|---|
| 1 | CALSPC | C*1 | Filename in ASCII original calibration file (32 characters) |
| 17 | CREDAT | C*1 | ASCII date calibration file was created (10 characters). The last character should be ignored if it is not a "1". If it is a "1", it indicates the data is after the year 2000. |
| 22 | CRETIM | C*1 | ASCII time calibration file was created (8 characters) |
| 26–32 | | | Reserved |
| 33 | EFFDAT | C*1 | ASCII date when spectrum for efficiency calibration was collected (10 characters) |
| 38 | EFFTIM | C*1 | ASCII time efficiency calibration file was created (8 characters) |
| 42–64 | | | Unused |

If long filenames are enabled, then the following 2 sequential records are defined as follows:

1+2 Filename in ASCII original calibration file

The energy channel pairs are stored as three records for peak centroid, three records for peak energy, and three records for peak FWHM. The three records each contain 32 four-byte real numbers for a total of 96 data points. The three sets of numbers are arranged in ascending energy order. The records are stored sequentially beginning with the record number in Word 22 of the first record. They are in the following order:

| <u>Word-22</u> | <u>Contents</u> |
|----------------|-------------------------------|
| +0 | Centroid values from 1 to 32 |
| +1 | Energy values from 1 to 32 |
| +2 | FWHM values from 1 to 32 |
| +3 | Centroid values from 33 to 64 |
| +4 | Energy values from 33 to 64 |
| +5 | FWHM values from 33 to 64 |
| +6 | Centroid values from 65 to 96 |
| +7 | Energy values from 65 to 96 |
| +8 | FWHM values from 65 to 96 |

The following is a description of an efficiency pairs record(s). The first such record is pointed to by Word 20 of the first record; the remaining records (if any) are positioned sequentially for the number of records in Word 30 of the first record.

| <u>Word Number</u> | <u>Local Name</u> | <u>Type</u> | <u>Use</u> |
|------------------------|-----------------------|-------------|--|
| 1 | EFF(1,1) | R*4 | Energy of first energy-efficiency pair in keV |
| 3 | EFF(2,1) | R*4 | Efficiency of first energy-efficiency pair |
| 5 | EFF(1,2) | R*4 | Energy of second energy-efficiency pair in keV |
| . | . | . | |
| . | . | . | |
| . | . | . | |
| 61 | EFF(1,16) | R*4 | Energy of 16th energy-efficiency pair in keV |
| 63 | EFF(2,16) | R*4 | Efficiency of 16th energy-efficiency pair |

4.3.2.2. Sodium Iodide Calibration

Record 1 (only valid if first analysis record, word 1 is 2 in the low byte)

| <u>Byte Number</u> | <u>Type</u> | <u>Use</u> |
|------------------------|-------------|--|
| 1 | AFIT | Above-knee efficiency calibration fit type (0 = None; 1 = Interp; 2 = Linear; 3 = Quad; 4 = Cubic) |
| 2 | DETNUM | Detector number |
| 3 | EFFPRS | # of efficiency pairs for interpolation |
| 4 | NCH | No. of channels in spectrum |
| 5–6 | ASIG | Error in fit above knee (2-sigma uncertainty) |
| 7–8 | | Error in FWHM at bottom channel |
| 9–10 | | Error in FWHM at top channel |
| 11–16 | | Energy calibration (offset, linear, quadratic parameters) |
| 17–22 | | FWHM calibration (offset, linear, quadratic parameters) |
| 23–30 | RE(1–4) | Cubic efficiency fit |
| 31 | ENGPRS | Number of energy triples |
| 32–33 | CALUNC | Calibration source uncertainty Record 2 (only valid if first analysis record, word 1 is 2 in the low byte) |

| <u>Byte Number</u> | <u>Type</u> | <u>Use</u> |
|------------------------|-------------|---|
| 1–128 | CALSPC | Filename in ASCII for original calibration file |
| 129–135 | CREDAT | ASCII calibration creation date |
| 136–140 | CRETIM | ASCII calibration creation time |

| | | |
|---------|--------|---|
| 141–146 | EFFDAT | ASCII date when spectrum for efficiency calibration was collected |
| 147–151 | EFFTIM | ASCII time when spectrum for efficiency calibration was collected |
| 152–184 | | Description line 1 |
| 184–216 | | Description line 2 |
| 217–512 | | Reserved |

4.4. Detector Description Files

A detector description file is created by the Detector program (A18) and contains the detector description input by the operator. The file extension should be **.DET**. These are not used in Windows programs. The contents of the file are:

Record 1

| Word Number | Use |
|----------------|-------------------------------------|
| 1 | Must be 1 |
| 2 | Must be 16 |
| 3–6 | Reserved |
| 7 | Detector description record pointer |
| 8–27 | Unused |
| 28 | Maximum record number ever used |
| 29 | Maximum record number in use |
| 30–50 | Reserved |
| 51 | Detector identification number |
| 52–64 | Reserved |

4.4.1. Detector Description Record

The detector description record contains the user-supplied description of the detector used to acquire the spectrum. The record is divided into two 64-character lines of text.

4.5. Old Geometry Correction Files

A geometry correction file is created by the program GEO (A30) and contains the correction information used in the analysis program (A30). The contents are:

Record 1

| <u>Word</u> <u>Number</u> | <u>Use</u> |
|------------------------------|--|
| 1 | Must be 1 |
| 2 | Must be 128 |
| 3–14 | Reserved |
| 15 | Geometry correction description record pointer |
| 16 | Geometry correction data record pointer |
| 17–27 | Reserved |
| 28 | Maximum record number ever used |
| 29 | Maximum record number in use |
| 30–64 | Reserved |

4.5.1. Geometry Correction Data Record

The geometry correction data record contains the data required to perform geometry correction. The record is as follows:

| <u>Word</u> <u>Number</u> | <u>Local</u> <u>Name</u> | <u>Type</u> | <u>Use</u> |
|------------------------------|-----------------------------|-------------|---------------------------------------|
| 1 | GEKEV(1) | R*4 | Energy of the first datum in keV |
| 3 | GEKEV(2) | R*4 | Energy of the second datum in keV |
| . | . | . | |
| . | . | . | |
| . | . | . | |
| 31 | GEKEV(16) | R*4 | Energy of the 16th datum in keV |
| 33 | GEVAL(1) | R*4 | Correction value for the first datum |
| 35 | GEVAL(2) | R*4 | Correction value for the second datum |
| . | . | . | |
| . | . | . | |
| . | . | . | |
| 63 | GEVAL(16) | R*4 | Correction value for the 16th datum |

4.5.2. Geometry Correction Description Record

The geometry correction description record contains information about the data on the geometry correction data record. Note that this record differs from the geometry record in the **.UFO** file (Section 4.15.5). The record is as follows:

| <u>Word Number</u> | <u>Local Name</u> | <u>Type</u> | <u>Use</u> |
|------------------------|-----------------------|-------------|--|
| 1 | NPAIRS | I*2 | Number of data pairs on data record |
| 2 | DTYPE | I*2 | 1 indicates geometry correction |
| 3 | GDETID | I*2 | Detector identification number |
| 4 | DATE1 | C*1 | Date geometry file was created (DD-MMM-YY*) (10 characters). The * character should be ignored if it is not a "1". If it is a "1", it indicates the data is after the year 2000. |
| 9 | TIME1 | C*1 | Time geometry file was created (HH:MM:SS) (8 characters) |
| 13 | FILESP | C*1 | Original geometry filename in ASCII (32 characters) |
| 29 | IDENT | C*1 | Detector description in ASCII (32 characters) |
| 49–64 | | | Reserved |

4.6. New Geometry Correction Files

The new geometry correction master file is written by GammaVision and contains the energy, Area1, Area2, ratio values, reference file, current file, and Detector ID. It has a file extension of **.GEO**. This file contains the following format, which is the same as the **.SOR** file format:

```

Line 1:      date and time
Line 1:      header
Lines 2–n:   Energy, Area1, Area2, Ratio values
blank line
Length:      (absorption only)
blank line   (absorption only)
Reference File:
Current File:
blank line
Detector ID:  (geometry only)

```

If long filenames are enabled, then the following 4 sequential records are defined as follows:

```

1+2          Reference UFO
3+4          Current UFO

```

4.7. Sample Description Files

A sample description file is written by the program SAMPLE (A18) and contains the sample description entered by the operator. It has a file extension of **.SMP**. This file is not used in GammaVision. This file contains the following:

Record 1

| Word | |
|---------------|-----------------------------------|
| <u>Number</u> | <u>Use</u> |
| 1 | Must be 1 |
| 2 | Must be 8 |
| 3–5 | Reserved |
| 6 | Sample description record pointer |
| 7–27 | Reserved |
| 28 | Maximum record number ever used |
| 29 | Maximum record number in use |
| 30–64 | Reserved |

4.7.1. Sample Description Record

This record is 128 ASCII characters split into two 64-character lines exactly as input by the operator. Unused characters are set ASCII spaces.

4.8. ROI File

The ROI (region-of-interest) file is created by the program ROIDEFIN (A18) and contains the definition of the regions of interest entered by the operator. The file contents are as follows:

4.8.1. Record 1

| Word | |
|---------------|---|
| <u>Number</u> | <u>Use</u> |
| 1 | Must be 1 |
| 2 | Must be 2048 |
| 3–20 | Reserved |
| 21 | Record number of the first of the two ROI records |
| 22–27 | Reserved |
| 28 | Maximum record number ever used |
| 29 | Maximum record number in use |
| 30–64 | Reserved |

4.8.2. ROI Data Record 1 and 2

There are two ROI data records, each record containing 64 short integer (2 byte) values. The first record contains the starting channel number for up to 64 ROIs. The second record, which always follows ROI Record 1, contains the number of channels for each of up to 64 ROIs. Undefined ROIs will have -1 as start channel number. Entries in Record 2 for undefined ROIs are undefined, but should be written as -1. The records are formatted as follows:

ROI Record 1

| Word Number | Type | Use |
|----------------|------|---------------------------------|
| 1 | I*2 | Start channel number for ROI 1 |
| 2 | I*2 | Start channel number for ROI 2 |
| 3 | I*2 | Start channel number for ROI 3 |
| . | . | . |
| . | . | . |
| . | . | . |
| 64 | I*2 | Start channel number for ROI 64 |

| Word Number | Type | Use |
|----------------|------|------------------------------|
| 1 | I*2 | Number of channels in ROI 1 |
| 2 | I*2 | Number of channels in ROI 2 |
| 3 | I*2 | Number of channels in ROI 3 |
| . | . | . |
| . | . | . |
| . | . | . |
| 64 | I*2 | Number of channels in ROI 64 |

4.9. Old Absorption Correction Files

The absorption correction file is written by the SOR program (A30) and contains the correction information used by the analysis programs. The contents of the file are as follows:

| Word Number | Use |
|----------------|--|
| 1 | Must be 1 |
| 3 | Must be 64 |
| 3-12 | Reserved |
| 13 | Absorption correction description record pointer |
| 14 | Absorption correction data record pointer |

| | |
|-------|---------------------------------|
| 15–27 | Reserved |
| 28 | Maximum record number ever used |
| 29 | Maximum record number in use |
| 30–64 | Reserved |

4.9.1. Absorption Correction Data Record

The absorption correction data record contains the data required to perform absorption correction. The record is formatted as follows:

| <u>Word</u> <u>Number</u> | <u>Local</u> <u>Name</u> | <u>Type</u> | <u>Use</u> |
|------------------------------|-----------------------------|-------------|---------------------------------------|
| 1 | SRKEV(1) | R*4 | Energy of the first datum in keV |
| 3 | SRKEV(2) | R*4 | Energy of the second datum in keV |
| . | . | . | . |
| . | . | . | . |
| . | . | . | . |
| 31 | SRKEV(16) | R*4 | Energy of the 16th datum in keV |
| 33 | SRVAL(1) | R*4 | Correction value for the first datum |
| 35 | SRVAL(2) | R*4 | Correction value for the second datum |
| . | . | . | . |
| . | . | . | . |
| . | . | . | . |
| 63 | SRVAL(16) | R*4 | Correction value for the 16th datum |

4.9.2. Absorption Correction Description Record

The absorption correction description record, indexed by Word 13 of the first record, contains information about the data contained on the absorption correction data record. See also the UFO Absorption record (Section 4.15.4). The record is formatted as follows:

| <u>Word</u> <u>Number</u> | <u>Local</u> <u>Name</u> | <u>Type</u> | <u>Use</u> |
|------------------------------|-----------------------------|-------------|--|
| 1 | NPAIRS | I*2 | Number of data pairs in data record |
| 2 | DTYPE | I*2 | 2 = internal, 4 = external |
| 3 | | C*2 | Reserved |
| 4 | DATE1 | C*10 | Date absorption file was created (DD-Mmm-YY); (stored as 10 ASCII characters). The * character should be ignored if it is not a "1". If it is a "1", it indicates the data is after the year 2000. |

| | | | |
|-------|--------|-----|--|
| 9 | TIME1 | C*1 | Time absorption file was created (HH:MM:SS); (stored as 8 ASCII characters) |
| 14 | | C*2 | Reserved |
| 16 | FILESP | C*1 | Original absorption (stored as 32 ASCII characters) |
| 32 | | C*2 | Reserved |
| 33 | IDENT | C*1 | Unit of measure as 30 ASCII characters (e.g., grams) |
| 49–64 | | | Reserved |

4.10. New Absorption Correction Files

The new absorption correction master file is written by GammaVision and contains the energy, area1, area2, ratio values, length, reference file, current file, and Detector ID. It has a file extension of **.SOR**. This file contains the following format, which is the same as the **.GEO** file format:

| | |
|-----------------|------------------------------------|
| Line 1: | date and time |
| Line 1: | header |
| Lines 2–n: | Energy, Area1, Area2, Ratio values |
| blank line | |
| Length: | (absorption only) |
| blank line | (absorption only) |
| Reference File: | |
| Current File: | |
| blank line | |
| Detector ID: | (geometry only) |

If long filenames are enabled, then the following 4 sequential records are defined as follows:

| | |
|-----|---------------|
| 1+2 | Reference UFO |
| 3+4 | Current UFO |

4.11. Real Format Spectrum Files

A real format spectrum file is created by the program CONVERT (A18) and the file extension should be **.SPC**. This file contains a spectrum and the associated information required to analyze the spectrum. The associated information is obtained from other files related to the MCA and segment that acquired the spectrum and are merged by CONVERT.

The contents of the real format spectrum files are:

4.11.1. Record 1

| <u>Word</u> <u>Number</u> | <u>Local</u> <u>Name</u> | <u>Type</u> | <u>Use</u> |
|------------------------------|-----------------------------|-------------|--|
| 1 | INFTYP | I*2 | Must be 1 |
| 2 | FILTYP | I*2 | Must be 5 |
| 3 | | | Contents flag Bit 0 = 1 for long filename Bit 1 = 1 for ZDT spectrum and ROI in file |
| 4 | | | Reserved |
| 5 | ACQIRP | I*2 | Acquisition information record pointer |
| 6 | SAMDRP | I*2 | Sample description record pointer |
| 7 | DETDRP | I*2 | Detector description record pointer |
| 8 | EBRDESC | I*2 | EBR description record pointer |
| 9 | ANARP1 | I*2 | First analysis parameters record pointer |
| 10 | ANARP2 | I*2 | Second analysis parameters record pointer |
| 11 | ANARP3 | I*2 | Third analysis parameters record pointer |
| 12 | ANARP4 | I*2 | Fourth analysis parameters record pointer |
| 13 | SRPDES | I*2 | Absorption correction description record pointer |
| 14 | IEQDESC | I*2 | IEQ description record pointer |
| 15 | GEODES | I*2 | Geometry correction description record pointer |
| 16 | MPCDESC | I*2 | MPC description record pointer |
| 17 | CALDES | I*2 | Calibration description record pointer |
| 18 | CALRP1 | I*2 | First calibration data record pointer |
| 19 | CALRP2 | I*2 | Second calibration data record pointer |
| 20 | EFFPRP | I*2 | Efficiency pairs record pointer (first record) |
| 21 | ROIRP1 | I*2 | Record number of the first of the two ROI records |
| 22 | | | Energy pairs record pointer |
| 23 | | | Number of energy pair records |
| 24 | | | Reserved |
| 25 | | | Disable deconvolution of unknown peaks |
| 26 | | | True = microcuries, false = becquerels |
| 27 | PERPTR | I*2 | Laboratory and operator name record pointer |
| 28 | MAXRCS | I*2 | Maximum record number ever used |
| 29 | LSTREC | I*2 | Maximum record number in use |
| 30 | EFFPNM | I*2 | Number of efficiency pairs records (See Word 20) |

| | | | |
|-------|--------|-----|---|
| 31 | SPCTRP | I*2 | Spectrum record pointer (pointer to first record) |
| 32 | SPCRCN | I*2 | Number of records in the spectrum |
| 33 | SPCCHN | I*2 | Number of channels in spectrum |
| 34 | ABSTCH | I*2 | Physical start channel for data |
| 35 | ACQTIM | R*4 | Date and time of acquisition start in DECDAY format |
| 37 | ACQTI8 | R*8 | Date and time as double precision DECDAY |
| 41 | SEQNUM | I*2 | Sequence number |
| 42 | MCANU | I*2 | MCA number as two ASCII characters (old) or Detector number as integer for systems with Connections |
| 43 | SEGNUM | I*2 | Segment number as two ASCII characters (old) or as integer value 1 for systems with Connections |
| 44 | MCADV | I*2 | MCA device type |
| 45 | CHNSRT | I*2 | Start channel number |
| 46 | RLTMDT | R*4 | Real Time in seconds |
| 48 | LVTMDT | R*4 | Live Time in seconds |
| 50 | | I*2 | Pointer to MGA or U235 or CZTU records |
| 51 | | I*2 | Pointer to FRAM records |
| 53 | | I*2 | Pointer to TRIFID records |
| 54 | | I*2 | Pointer to NaI records |
| 55–62 | | | Reserved |
| 63–64 | RRSFCT | R*4 | Total random summing factor |

4.11.2. Acquisition Information Record

The acquisition information record contains character data that is written by CONVERT from information in the **.CHN** file.

| Byte Number | Type | Use |
|----------------|------|--|
| 1 | C*1 | Default spectrum file name (stored as 16 ASCII characters) |
| 17 | C*1 | Date in the form DD-MMM-YY* (stored as 12 ASCII characters). The * character should be ignored if it is not a "1". If it is a "1", it indicates the data is after the year 2000. |
| 29 | C*1 | Time in the form HH:MM:SS (stored as 10 ASCII characters) |
| 39 | C*1 | Live Time rounded to nearest second (stored as 10 ASCII characters) |
| 49 | C*1 | Real Time rounded to nearest second (stored as 10 ASCII characters) |
| 59–90 | | Reserved |

| | | |
|-----|-----|---|
| 91 | C*1 | Start date of sample collection (10 ASCII characters) |
| 103 | C*1 | Start time of sample collection (8 ASCII characters) |
| 111 | C*1 | Stop date of sample collection (10 ASCII characters) |
| 121 | C*1 | Stop time of sample collection (8 ASCII characters) |

4.11.3. Sample Description Record

This record is discussed in Section 4.7.

4.11.4. Detector Description Record

This record is described in Section 4.4.

4.11.5. First, Second, and Third Analysis Parameter

These records are described in Section 4.2.

4.11.6. Absorption Correction Description Record

| <u>Word</u> <u>Number</u> | <u>Local</u> <u>Name</u> | <u>Type</u> | <u>Use</u> |
|------------------------------|-----------------------------|-------------|--|
| 1 | SORPRS | I*2 | Number of pairs (between 1 and 100) |
| 2 | SORINT | I*2 | 1 if geometry record, 2 if internal, 4 if external absorption |
| 3 | SORGEO | I*2 | 1 if cylindrical geometry (internal only), 2 if Marinelli geometry (internal only) |
| 4-5 | SORDEN | R*4 | Absorber density |
| 6-7 | SORLNG | R*4 | Absorber length |
| 8-12 | SORDAT | C*10 | Date absorption file created |
| 13-16 | SORTIM | C*10 | Time absorption file created |
| 17-32 | SORNAM | C*32 | Absorber name |
| 33-48 | SREFNAM | C*32 | Reference UFO name |
| 49-64 | SCURNAM | C*32 | Current UFO name |

If long filenames are enabled, then the following 4 sequential records are defined as follows:

| | |
|-----|--------------------|
| 1+2 | Reference UFO name |
| 3+4 | Current UFO name |

4.11.7. Absorption Correction Data Record

| Word Number | Local Name | Type | Use |
|----------------|----------------|------|---------------------------|
| 1–16 | SORKEV (1–8) | R*4 | Energy |
| 17–32 | SORAREA1 (1–8) | R*4 | Reference area |
| 33–48 | SORAREA2 (1–8) | R*4 | Current area |
| 49–64 | SORVAL (1–8) | R*4 | Linear attenuation values |

If fewer than 8 pairs are used, the array values are set to zero. There are as many data records as necessary, up to 64 energies.

4.11.8. Geometry Correction Description Record

| Word Number | Local Name | Type | Use |
|----------------|---------------|------|---|
| 1 | GEOPRS | I*2 | Number of pairs (between 1 and 100) |
| 2 | GEOINT | I*2 | 1 if geometry record, 2 if internal, 4 if external absorption |
| 3 | unused | I*2 | |
| 4–5 | unused | R*4 | |
| 6–7 | unused | R*4 | |
| 8–12 | GEODAT | C*10 | Date geometry file created |
| 13–16 | GEOTIM | C*10 | Time geometry file created |
| 17–32 | unused | C*32 | |
| 33–48 | GREFNAM | C*32 | Reference UFO name |
| 49–64 | GCURNAM | C*32 | Current UFO name |

If long filenames are enabled, then the following 4 sequential records are defined as follows:

| | |
|-----|--------------------|
| 1+2 | Reference UFO name |
| 3+4 | Current UFO name |

4.11.9. Geometry Correction Data Record

| Word Number | Local Name | Type | Use |
|----------------|----------------|------|----------------------------|
| 1–16 | GEOKEV (1–8) | R*4 | Energy |
| 17–32 | GEOAREA1 (1–8) | R*4 | Reference peak area |
| 33–48 | GEOAREA2 (1–8) | R*4 | Current peak area |
| 49–64 | GEOVAL (1–8) | R*4 | Geometry correction ratios |

If fewer than 8 pairs are used, the array values are set to zero. Up to 64 energies can be stored. The unused values are set to zero.

4.11.10. Calibration Description Record, Calibration Data Records, and Efficiency Pairs Records

These records are described in Section 4.3.

4.11.11. Region-of-Interest Records

This format follows the format described in Section 3.2. The ROI records will be three consecutive records, with the first record pointed to by word 21 of record 1 (FORTRAN name = ROIRP1). The end of the list of regions is given by a negative start channel.

Record 1

| <u>Word</u> <u>Number</u> | <u>Use</u> |
|------------------------------|-----------------------------------|
| 0 | Must be -2 |
| 1 | Start channel of first ROI |
| 2 | Stop channel of first ROI |
| 3 | Start channel of second ROI |
| 4 | Stop channel of second ROI |
| . | . |
| . | . |
| . | . |
| 61 | Start channel of thirty-first ROI |
| 62 | Stop channel of thirty-first ROI |
| 63 | Reserved |

A negative start channel terminates the list of regions.

Record 2

| <u>Word</u> <u>Number</u> | <u>Use</u> |
|------------------------------|------------------------------------|
| 0 | Start channel of thirty-second ROI |
| 1 | Stop channel of thirty-second ROI |
| . | . |
| . | . |
| . | . |

- 62 Start channel of 63rd ROI
63 Stop channel of 63rd ROI

Record 3

This contains the energy calibration for the MCB or buffer at the time the ROIs were defined.

| Word Number | Type | Use |
|----------------|------|-----------------------------------|
| 0 | R*4 | Energy calibration zero intercept |
| 3 | R*4 | Energy calibration slope |
| 7 | R*4 | Energy calibration quadratic term |
| 11–63 | | Reserved |

4.11.12. Hardware Parameters Records

These records contain the hardware parameters for the MCB when the data were collected. The appropriate values are entered into the record when the spectrum file is saved. Any value that is correctly inserted must have the corresponding bit set in the validity flag words at the beginning of the first record. All unused bytes are set to binary 0. Any program reading this record must check the validity bit for the specific value before using the value in the record.

These two sequential records are pointed to by word 12 of record 1 in any inform file (FORTRAN name = ANARP4).

Record 1

| Word Number | Type | Use |
|----------------|------|--|
| 1–4 | | Validity flag, each bit corresponds to a single entry. 1 = valid contents, 0 = unused for this record. Bits are counted from the right, starting with 1. The bit number corresponds to the word number; i.e., if preset real time (word 5) is valid, then bit 5 and 6 are 1. |
| 5 | R*4 | Preset real time in seconds |
| 7 | R*4 | Preset live time in seconds |
| 9 | I*4 | Preset counts |
| 11 | I*4 | Preset integral |
| 13 | R*4 | Preset uncertainty in percent |
| 15 | | Overflow preset, 0 = off, 1 = on |
| 16 | I*2 | Start channel for uncertainty preset region |
| 17 | I*2 | Stop channel for uncertainty preset region |
| 18 | I*2 | ADC conversion gain in channels |

| | | |
|----|-----|--|
| 19 | I*2 | ADC offset in channels |
| 20 | I*2 | Lower level discriminator in channels |
| 21 | I*2 | Upper level discriminator in channels |
| 22 | I*2 | Input gate: 0 = off, 1 = coincidence, 2 = anti |
| 23 | I*2 | Amplifier coarse gain |
| 24 | R*4 | Amplifier fine gain as the value of the multiplier seen on the MAESTRO display |
| 26 | R*4 | Amplifier fine offset in fractional channels |
| 28 | I*2 | Gain stabilizer adjustment amplifier gain setting as the value from the MCB |
| 29 | R*4 | Start channel for gain stabilizer region |
| 31 | R*4 | Stop channel for gain stabilizer region |
| 33 | | Gain stabilizer mode; 0 = Gauss, 1 = peak |
| 34 | | Gain stabilizer on/off; 1 = on, 0 = off |
| 35 | I*2 | Zero stabilizer adjustment setting as the value from the MCB |
| 36 | R*4 | Start channel for zero-stabilizer region |
| 38 | R*4 | Stop channel For zero-stabilizer region |
| 40 | | Zero stabilizer on/off; 1 = on, 0 = off |
| 41 | R*4 | Shaping time constant in microseconds (as reported) |
| 43 | I*2 | Preamplifier type; 0 = resister, 1 = TRP |
| 44 | I*2 | PZ valid; 0 = no, 1 = yes |
| 45 | I*2 | PZ value as reported by MCB |
| 46 | I*2 | High voltage value in volts |
| 47 | R*4 | Rise time (DSPEC [®]) as reported on MAESTRO display |
| 49 | R*4 | Width (DSPEC) as reported on MAESTRO display |
| 51 | R*4 | Cusp (DSPEC) as reported on MAESTRO display |
| 53 | R*4 | Tilt (DSPEC) as reported on MAESTRO display |
| 55 | I*2 | Baseline type; 1 = auto, 2 = fast, 3 = slow |
| 56 | R*4 | Baseline value in microseconds |
| 58 | I*2 | Amplifier mode, 1 = germanium, 2 = sodium iodide |
| 59 | I*2 | Acquisition mode, 1 = PHA, 2 = MCS |
| 60 | | MCB serial number, as ASCII text, 8 characters |
| 64 | I*2 | Number of hardware records in file |

Record 2

| Word | | |
|---------------|-------------|--|
| <u>Number</u> | <u>Type</u> | <u>Use</u> |
| 1–2 | I*4 | Validity flags for words 1–32, each bit corresponds to a single entry. 1 = valid contents, 0 = unused for this record |
| 3–4 | I*4 | Validity flags for words 33–64 |

| | | |
|-------|-------|--|
| 5 | I*2 | Amplifier input polarity, +1 = positive, -1 = negative |
| 6–7 | R*4 | Thermistor value, in ohms |
| 8 | I*2 | Pileup rejector, 0 = off, 1 = on |
| 9–10 | I*4 | Pileup rejector width in ns |
| 11 | I*2 | OCTÊTE-PC current in nA |
| 12 | I*2 | OCTÊTE-PC vacuum in mT |
| 13–20 | C*8 | Firmware revision in ASCII, 16 characters |
| 21 | I*2 | High Voltage enabled: 1 = Yes, 0 = No |
| 22 | I*2 | Shaping time constant index |
| 23 | I*2 | Battery Stat: 0 = Ext, 1 = Battery 1, 2 = Battery 2 |
| 24–25 | R*4 | Battery 1 voltage |
| 26–27 | R*4 | Battery 2 voltage |
| 28 | I*2 | Desired high voltage |
| 29 | I*2 | HV shutdown mode: 0 = TTL, 1 = ORTEC, 2 = Off |
| 30 | I*2 | ADC type: 0 = CI34, 1 = CI36, 2 = ORTEC, 3 = Silena |
| 31 | I*2 | Maximum vacuum value |
| 32 | I*2 | Automatic threshold flag (SBS-60): 0 = Off, 1 = On |
| 33–38 | float | MDA preset coefficients |
| 39–40 | float | MDA preset value |
| 41–42 | float | User's MDA preset value |
| 43 | I*2 | Low MDA preset ROI limit |
| 44 | I*2 | High MDA preset ROI limit |
| 45 | short | MDA preset units: 0 = Bq, 1 = μ Ci, -1 = None |
| 46–49 | char | MDA preset nuclide name in ASCII |
| 50 | I*2 | ZDT-enabled flag: 1 = On, 0 = Off |
| 51 | I*2 | ZDT refresh rate in microsecs |
| 52 | I*2 | ZDT view spectrum (1 = Normal, 2 = Corrected) |
| 53–54 | R*8 | Total count rate in CPS |
| 55 | I*2 | MiniMCA-166: analog threshold value |
| 56 | I*2 | MiniMCA-166: input routing 0 = through internal amplifier 1 = bypass internal amplifier with 0 to +3 V 2 = bypass internal amplifier with 0 to -3 V |
| 57 | I*2 | Firmware revision code (MiniMCA-166) |
| 58 | I*2 | Hardware revision code (MiniMCA-166) |
| 59 | I*2 | Time value low for PZ (MiniMCA-166) |
| 60 | I*2 | Time value high for PZ (MiniMCA-166) |
| 61 | I*2 | Slow discriminator value (MiniMCA-166) |
| 62 | I*2 | Fast discriminator value (MiniMCA-166) |

| | | |
|----|-----|--|
| 63 | I*2 | Power switches (MiniMCA-166) |
| 64 | I*2 | ZDT mode: 0 = Off, 1 = Normal-Corrected (LTC-ZDT), 2 = Error-Uncertainty (ERR-ZDT) |

Record 3

| <u>Word</u> <u>Number</u> | <u>Local</u> <u>Name</u> | <u>Type</u> | <u>Use</u> |
|------------------------------|-----------------------------|-------------|--|
| 1–2 | dwFlg1 | I*4 | Validity flags for words 1–32 |
| 3–4 | dwFlg2 | I*4 | Validity flags for words 33–64 |
| 5 | wP12I | I*2 | +12-V current in mA (MiniMCA-166) |
| 6 | wM12I | I*2 | –12-V current in mA (MiniMCA-166) |
| 7 | wP24I | I*2 | +24-V current in mA (MiniMCA-166) |
| 8 | wM24I | I*2 | –24-V current in mA (MiniMCA-166) |
| 9 | wBatI | I*2 | Battery current in mA (MiniMCA-166) |
| 10 | wHvI | I*2 | HV current in mA (MiniMCA-166) |
| 11 | wChrI | I*2 | Charger current in mA (MiniMCA-166) |
| 12 | wColl | I*2 | Collection time (SBS-60) |
| 13 | wCollFast | I*2 | Collection start on fast signal: 0 = No, 1 = Yes |
| 14 | wCollSlow | I*2 | Collection start on slow signal: 0 = No, 1 = Yes |
| 15 | wHVRise | I*2 | HV rise time in seconds |
| 16 | wFastThr | I*2 | Fast threshold value (SBS-60) |
| 17 | wSyncT | I*2 | Synchronize thresholds flag: : 0 = No, 1 = Yes |
| 18 | wDifTC | I*2 | Differentiation time constant in μ s |
| 19 | wIntTC | I*2 | Integration time constant in μ s |
| 20 | wInhbt | I*2 | Inhibit signal polarity (SBS-60): 0 = Off, 1 = High, 2 = Low |
| 21 | | | Number of valid nuclides in MDA preset |
| 22 | | | High voltage polarity ± 1 |
| 22–30 | | | Text associated with shutdown modes |
| 31–32 | | | HV current in mA (good for all HVPS) |
| 33–40 | | | Text associated ADC type |
| 41–48 | | | Text associated with Gate |
| 49–56 | | | Text associated with PZ mode |
| 57 | | | Number of strings with szViewNames |
| 58 | | | HV shutdown status () = OK, 1 = Shutdown) |
| 59–60 | | | Gain stabilizer adjustment in percent (–100 to +100) |
| 61–62 | | | Zero stabilizer adjustment in percent (–100 to +100) |
| 63 | | | Number of status string pairs following View strings |
| 64 | | | Unused |

Record 4

| Word | |
|---------------|-----------------------------------|
| <u>Number</u> | <u>Use</u> |
| 1–2 | Validity flags for 1–32 |
| 3–4 | Validity flats for 33–64 |
| 5 | Start delay in seconds (DART) |
| 6 | Conserve delay in seconds (DART) |
| 7 | Off delay in seconds (DART) |
| 8 | Power mode (0 = On, 1 = Conserve) |
| 9–32 | Unused |
| 33–64 | Data view name strings |

Multi-Nuclide MDA Preset Storage

Up to 20 nuclides according to Word 21 in record 3, which corresponds to records 5–9.

| Word | | |
|---------------|-------------|-------------------------|
| <u>Number</u> | <u>Type</u> | <u>Use</u> |
| 1–2 | R*4 | MDA user preset (in Bq) |
| 3–4 | R*4 | Eff * Yield |
| 5–6 | R*4 | MDA hardware preset |
| 7–8 | R*4 | MDA energy |
| 9 | I*2 | Low ROI limit |
| 10 | I*2 | High ROI limit |
| 11 | | Reserved |
| 12–16 | C*5 | MDA nuclide 1 name |
| 17–18 | | MDA user preset (in Bq) |
| 19–20 | | Eff * Yield |
| 21–22 | | MDA hardware preset |
| 23–24 | | MDA energy |
| 25 | | Low ROI limit |
| 26 | | High ROI limit |
| 27 | | Reserved |
| 28–32 | | MDA nuclide 2 name |
| 33–34 | | MDA user preset (in Bq) |
| 35–36 | | Eff * Yield |
| 37–38 | | MDA hardware preset |
| 39–40 | | MDA energy |
| 41 | | Low ROI limit |
| 42 | | High ROI limit |
| 43 | | Reserved |

| | |
|-------|-------------------------|
| 44–48 | MDA nuclide 3 name |
| 49–50 | MDA user preset (in Bq) |
| 51–52 | Eff * Yield |
| 53–54 | MDA hardware preset |
| 55–56 | MDA energy |
| 57 | Low ROI limit |
| 58 | High ROI limit |
| 59 | Reserved |
| 60–64 | MDA nuclide 4 name |

This is repeated in blocks of 4 nuclides for the total needed. Records 10 and up contain status records for the hardware. The number of records is Word 63 of record 3.

4.11.13. Personality Record

This record contains the laboratory name and the operator name for this spectrum. These are set in the GammaVision dialogs.

| Word | Local | | |
|---------------|-------------|-------------|-----------------|
| <u>Number</u> | <u>Name</u> | <u>Type</u> | <u>Use</u> |
| 1 | LABNAM | C*64 | Laboratory name |
| 33 | OPRNAM | C*64 | Operator name |

4.11.14. Table Description Record

This set of records is pointed to by the TABREC pointer of record 1. It contains the tables for the IEQ, EBAR, and MPC corrections.

| Word | Local | | |
|---------------|-------------|-------------|--------------------------------------|
| <u>Number</u> | <u>Name</u> | <u>Type</u> | <u>Use</u> |
| 1 | ITBLTYP | I*2 | Data type: 3 = EBR, 4 = IEQ, 5 = MPC |
| 2 | NTBLVAL | I*2 | Number of data pairs |
| 3–7 | TBLDAT | C*10 | Date geometry file created |
| 8–11 | TBLTIM | C*8 | Time geometry file created |
| 12–18 | MPCUNIT | C*14 | MPC data units |
| 19–34 | TABNAM | C*32 | Table filename |

If long filenames are enabled, then the following 2 sequential records are defined as follows:

| | |
|-----|----------------|
| 1+2 | Table filename |
|-----|----------------|

4.11.15. Table Data Record

This record follows the table description record. Any unused values are set to -1.

| Word | Local | | |
|---------------|---------------|-------------|--|
| <u>Number</u> | <u>Name</u> | <u>Type</u> | <u>Use</u> |
| 1–40 | ISONAM (1–10) | C*8 | Array of isotope names |
| 41–61 | TABVAL (1–10) | R*4 | Array of average energy or IEQ or MPC values |

Table Data Records 2–N follow if necessary. Up to 40 isotopes can be stored. The unused values are set to zero.

4.11.16. Spectrum Records

These records are the spectrum data stored as REAL*4 numbers beginning with the channel number given and going through the number of channels in the file. They are stored as 64-word records, which gives 32 data channels per record. They are stored sequentially, beginning with the record pointer given.

4.12. Integer Format Spectrum Files

An integer format spectrum file is created by the GammaVision and the file extension should be **.SPC**. This file contains a spectrum and the associated information required to analyze the spectrum. The contents of the integer format spectrum files are as follows.

Record 1

| Word | Local | | |
|---------------|-------------|-------------|---|
| <u>Number</u> | <u>Name</u> | <u>Type</u> | <u>Use</u> |
| 1 | INFTYP | I*2 | Must be 1 |
| 2 | FILTYP | I*2 | Must be 1 |
| 3 | | | Contents flag Bit 0 = 1 for long filename Bit 1 = 1 for ZDT spectrum and ROI in file. |
| 4 | | | Reserved |
| 5 | ACQIRP | I*2 | Acquisition information record pointer |
| 6 | SAMDRP | I*2 | Sample description record pointer |
| 7 | DETDRP | I*2 | Detector description record pointer |
| 8 | EBRDESC | I*2 | EBAR description record |
| 9 | ANARP1 | I*2 | First analysis parameters record pointer |
| 10 | ANARP2 | I*2 | Second analysis parameters record pointer |
| 11 | ANARP3 | I*2 | Third analysis parameters record pointer |

| | | | |
|----|---------|-----|---|
| 12 | ANARP4 | I*2 | Fourth analysis parameters record pointer |
| 13 | SRPDES | I*2 | Absorption correction description record pointer |
| 14 | IEQDESC | I*2 | IEQ description record pointer |
| 15 | GEODES | I*2 | Geometry correction description record pointer |
| 16 | MPCDESC | I*2 | MPC description record pointer |
| 17 | CALDES | I*2 | Calibration description record pointer |
| 18 | CALRP1 | I*2 | First calibration data record pointer |
| 19 | CALRP2 | I*2 | Second calibration data record pointer |
| 20 | EFFPRP | I*2 | Efficiency pairs record pointer (first record) |
| 21 | ROIRP1 | I*2 | Record number of the first of the two ROI records |
| 22 | | | Energy pairs record pointer |
| 23 | | | Number of energy pair records |
| 24 | | | Reserved |
| 25 | | | Disable deconvolution of unknown peaks |
| 26 | | | True = microcuries, false = becquerels |
| 27 | PERPTR | I*2 | Laboratory and operator name record pointer |
| 28 | MAXRCS | I*2 | Maximum record number ever used |
| 29 | LSTREC | I*2 | Maximum record number in use |
| 30 | EFFPNM | I*2 | Number of efficiency pairs records (See Word 20) |
| 31 | SPCTRP | I*2 | Spectrum record pointer (pointer to first record) |
| 32 | SPCRCN | I*2 | Number of records in the spectrum |
| 33 | SPCCHN | I*2 | Number of channels in spectrum |
| 34 | ABSTCH | I*2 | Physical start channel for data |
| 35 | ACQTIM | R*4 | Date and time of acquisition start in DECDAY format |
| 37 | ACQTI8 | R*8 | Date and time as double precision DECDAY |
| 41 | SEQNUM | I*2 | Sequence number |
| 42 | MCANU | I*2 | MCA number as two ASCII characters (old) or Detector number as integer for systems with Connections |
| 43 | SEGNUM | I*2 | Segment number as two ASCII characters (old) or as integer value 1 for systems with Connections |
| 44 | MCADVT | I*2 | MCA device type |
| 45 | CHNSRT | I*2 | Start channel number |
| 46 | RLTMDT | R*4 | Real time in seconds |
| 48 | LVTMDT | R*4 | Live time in seconds |

50
51–64

Reserved
Reserved

4.12.1. Other Records

The remaining records are described in the real format spectrum file section with the exception that the spectrum is stored as 4-byte integers.

4.12.2. Spectrum Records

These records are the spectrum data stored as INTEGER*4 numbers beginning with the channel number given and going through the number of channels in the file. They are stored as 64-word records, which gives 32 data channels per record. They are stored sequentially, beginning with the record pointer given.

4.13. Net Spectrum Files

A net spectrum file should have the extension of **.SPC**. This file contains the difference between two spectra with the analysis parameters and other associated data taken from the foreground spectrum. The contents of the net spectrum file are:

Record 1

| <u>Word Number</u> | <u>Local Name</u> | <u>Type</u> | <u>Use</u> |
|------------------------|-----------------------|-------------|--|
| 1 | INFTYP | I*2 | Must be 1 |
| 2 | FILTYP | I*2 | Must be 7 (real) or 3 (integer) |
| 3 | | | Contents Flag Bit 0 = 1 for long filename Bit 1 = 1 for ZDT spectrum and ROI in file |
| 4 | | | Reserved |
| 5 | ACQIRP | I*2 | Acquisition information record pointer |
| 6 | SAMDRP | I*2 | Sample description record pointer |
| 7 | DETDRP | I*2 | Detector description record pointer |
| 8 | EBRDESC | I*2 | EBAR description record |
| 9 | ANARP1 | I*2 | First analysis parameters record pointer |
| 10 | ANARP2 | I*2 | Second analysis parameters record pointer |
| 11 | ANARP3 | I*2 | Third analysis parameters record pointer |
| 12 | ANARP4 | I*2 | Fourth analysis parameter record pointer |
| 13 | SRPDES | I*2 | Absorption correction description record pointer |

| | | | |
|-------|---------|-----|---|
| 14 | IEQDESC | I*2 | IEQ description data record pointer |
| 15 | GEODES | I*2 | Geometry correction description record pointer |
| 16 | MPCDESC | I*2 | MPC description record pointer |
| 17 | CALDES | I*2 | Calibration description record pointer |
| 18 | CALRP1 | I*2 | First calibration data record pointer |
| 19 | CALRP2 | I*2 | Second calibration data record pointer |
| 20 | EFFPRP | I*2 | Efficiency pairs record pointer (first record) |
| 21 | ROIRPI | I*2 | Record number of the first of the two ROI records |
| 22 | PERPTR | I*2 | Personality record pointer |
| 27 | CALRP2 | I*2 | Second calibration data record pointer |
| 28 | MAXRCS | I*2 | Maximum record number ever used |
| 29 | LSTREC | I*2 | Maximum record number in use |
| 30 | EFFPNM | I*2 | Number of efficiency pairs records (see Word 20) |
| 31 | SPCTRP | I*2 | Spectrum record pointer (pointer to first record) |
| 32 | SPCRCN | I*2 | Number of records in the spectrum |
| 33 | SPCCHN | I*2 | Number of channels in spectrum |
| 35 | ACQTIM | R*4 | Date and time of acquisition start in DECDAY format |
| 37 | ACQTI8 | R*8 | Date and time as double precision DECDAY |
| 41 | SEQNUM | I*2 | Sequence number |
| 42 | SEQNAM | C*1 | Sequence name of 4 ASCII characters |
| 44 | MCADVT | I*2 | MCA device type |
| 45 | CHNSRT | I*2 | Start channel number |
| 46 | RLTMDT | R*4 | Real time in seconds |
| 48 | LVTMDT | R*4 | Live time seconds |
| 50–64 | | | Reserved |

The remaining records are described in the real format spectrum file section.

4.14. AlphaVision Spectrum Files

These files are written by AlphaVision (A36-B32). They can be read by other programs (e.g., MAESTRO and WINPLOTS), but only the analysis parameters are used by AlphaVision.

The first record is:

| <u>Word</u> | <u>Use</u> |
|---------------|--|
| <u>Number</u> | |
| 1 | Must be 1 |
| 2 | Must be 1 |
| 3 | Contents flag; on is true |
| | <u>Bit</u> <u>Description</u> |
| | Bit 1 Long file names |
| 4 | Reserved |
| 5 | Acquisition information record pointer |
| 6 | Sample description record pointer |
| 7 | Detector description record pointer |
| 8 | First table record pointer |
| 9 | First analysis parameters record pointer |
| 10 | Second analysis parameters record pointer |
| 11 | Third analysis parameters record pointer |
| 12 | Fourth analysis parameters record pointer |
| 13 | Absorption correction description record pointer |
| 14 | Absorption correction data record pointer |
| 15 | Geometry correction description record pointer |
| 16 | Geometry correction data record pointer |
| 17 | Calibration description record pointer |
| 18 | First calibration data record pointer |
| 19 | Second calibration data record pointer |
| 20 | Efficiency pairs record pointer (first record) |
| 21 | ROI record pointer (first of the three records) |
| 22 | Energy pairs record pointer (first record) |
| 23 | Number of energy-pairs records |
| 24 | Reserved |
| 25 | Reserved |
| 26 | True = microcuries, False = becquerels |
| 27 | Pointer to personality record (laboratory name) |
| 28 | Maximum record allocated (unused for UFO) |
| 29 | Maximum record in use (or allocated for UFO) |
| 30 | Reserved |
| 31 | Spectrum record pointer |
| 32 | Number of records in spectrum |
| 33 | Number of channels |
| 34 | Start channel |
| 35–36 | Acquisition start in DECDAY format |

| | |
|-------|-------------------------------------|
| 37–40 | Acquisition start in DECDAY8 format |
| 41 | Sequence number |
| 42 | Detector ID |
| 43 | Segment number (always 1) |
| 44 | MCA device type |
| 45 | Start channel number |
| 46 | Real time in seconds |
| 48 | Live time in seconds |
| 50–64 | Reserved |

The remaining records are the same as the germanium spectrum files (see Section 4.12).

4.15. Germanium Unformatted Output Files

An unformatted output (.UFO) file is created by the programs AN1, FMPC (B30), MAN1, WAN1, GAM32 (A66), MFMP (A34), and all other Windows germanium-analysis programs. Older .UFO files are used by PBC and PEAKPLOT, but .UFO files created by GammaVision cannot be used with older programs. It should have a file extension of .UFO. A .UFO file contains the results, both intermediate and final, of the germanium spectral analysis. Because of their size, the user should delete .UFO files that are no longer needed.

The .UFO file has changed in GammaVision 5.2, that is, the “old” .UFO file was written before version 5.2 and the “new” .UFO file is used by version 5.2 and later. To determine if a .UFO file is “old” or “new,” look at Words 50 and 51 of Record 2. These will be **2** for new .UFO files and **0** (zero) for old files. Other values (i.e., not 2 or zero) should be interpreted as zero (old).

4.15.1. Records 1 and 2

4.15.1.1. Old UFO Records 1 and 2

The contents of the old .UFO file are:

Record 1

| <u>Word</u> <u>Number</u> | <u>Local</u> <u>Name</u> | <u>Use</u> |
|------------------------------|-----------------------------|--|
| 1 | | Must be 1 |
| 2 | | Must be 1024 |
| 3 | | Special application record pointer – 1 |
| 4 | | Special application record pointer – 2 |
| 5 | | Acquisition information record pointer |
| 6 | | Sample description record pointer |

| | | |
|-------|--------|--|
| 7 | | Detector description record pointer |
| 8 | | Special application record pointer |
| 9 | | First analysis parameters record pointer |
| 10 | | Second analysis parameters record pointer |
| 11 | | Third analysis parameters record pointer |
| 12 | | Fourth analysis parameters record pointer |
| 13 | | Absorption correction description record pointer |
| 14 | | Absorption correction data record pointer |
| 15 | | Geometry correction description record pointer |
| 16 | | Geometry correction data record pointer |
| 17 | | Calibration description record pointer |
| 18 | | First calibration data record pointer |
| 19 | | Second calibration data record pointer |
| 20 | | Efficiency pairs record pointer (first record) |
| 21 | | Record number of the first of the ROI records |
| 22 | | Energy pairs record pointer |
| 23 | | Number of energy pair records |
| 24 | | Reserved |
| 25 | | Disable deconvolution of unknown peaks |
| 26 | | True = microcuries, False = becquerels |
| 27 | PERPTR | Laboratory and operator name record pointer |
| 28 | | Maximum record number ever used |
| 29 | | Maximum record number in use |
| 30 | | Number of efficiency pairs records (see Word 20) |
| 31–64 | | Reserved |

Record 2

The second record of the old **.UFO** file contains pointers to other records used in the **.UFO** file and data used by the various analysis programs used to analyze the data.

| Word Number | Local Name | Use |
|----------------|---------------|---------------------------------------|
| 1 | I*2 | GEN record pointer |
| 2 | I*2 | CSI 1 record pointer |
| 3 | I*2 | CSI 2 record pointer |
| 4 | I*2 | Library 1 peak records start pointer |
| 5 | I*2 | Number of Library 1 peaks used |
| 6 | I*2 | Number of Library 1 peaks allocated |
| 7 | I*2 | Library 2 peak records start pointer |
| 8 | I*2 | Number of Library 2 peak records used |

| | | |
|----|-----|--|
| 9 | I*2 | Number of Library 2 peak records allocated |
| 10 | I*2 | Library 3 peak records start pointer |
| 11 | I*2 | Number of Library 3 peak records used |
| 12 | I*2 | Number of Library 3 peak records allocated |
| 13 | I*2 | Unknown peak records start pointer |
| 14 | I*2 | Number of unknown peak records used |
| 15 | I*2 | Number of unknown peak records allocated |
| 16 | I*2 | Library 1 nuclide records start pointer |
| 17 | I*2 | Number of Library 1 nuclide records used |
| 18 | I*2 | Number of Library 1 nuclide records allocated |
| 19 | I*2 | Library 2 nuclide records start pointer |
| 20 | I*2 | Number of Library 2 nuclide records used |
| 21 | I*2 | Number of Library 2 nuclide records allocated |
| 22 | I*2 | Library 3 nuclide records start pointer |
| 23 | I*2 | Number of Library 3 nuclide records used |
| 24 | I*2 | Number of Library 3 nuclide records allocated |
| 25 | | Reserved |
| 28 | I*2 | CSI 3 record pointer |
| 29 | I*2 | CSI 4 record pointer |
| 30 | I*2 | Lowest in-range peak pointer for Lib-1 and this analysis |
| 31 | I*2 | Lowest in-range peak pointer for Lib-2 and this analysis |
| 32 | I*2 | Lowest in-range peak pointer for Lib-3 and this analysis |
| 33 | R*4 | NUPISA |
| 35 | R*4 | EBRDEC |
| 37 | R*4 | IEQVAL |
| 39 | R*4 | IEQDEC |
| 41 | I*2 | NAC shadow for library 1 |
| 43 | I*2 | NAC shadow for library 2 |
| 45 | I*2 | NAC shadow for library 3 |
| 49 | | Pointer to start record 3 |
| 50 | | Reserved |
| 59 | R*4 | EBRVAL |
| 61 | R*4 | Total activity |
| 63 | R*4 | Etime |

Other Records

The records not described below are described in the real format spectrum file discussion, Section 4.11.

4.15.1.2. New UFO Records 1 and 2

The contents of the new .UFO file are:

Record 1

| <u>Word Number</u> | <u>Local Name</u> | <u>Use</u> |
|------------------------|-----------------------|---|
| 1 | | Must be 1 |
| 2 | | Must be 1024 |
| 3 | | Contents flag, 1 = true |
| | | <u>Bit</u> <u>Description</u> |
| | | Bit 0 Long filename |
| | | Bit 1 ZDT |
| 4 | | Reserved |
| 5 | | Acquisition information record pointer |
| 6 | | Sample description record pointer |
| 7 | | Detector description record pointer |
| 8 | | First table record pointer |
| 9 | | First analysis parameters record pointer |
| 10 | | Second analysis parameters record pointer |
| 11 | | Third analysis parameters record pointer |
| 12 | | Fourth analysis parameters record pointer |
| 13 | | Absorption correction description record pointer |
| 14 | | Absorption correction data record pointer |
| 15 | | Geometry correction description record pointer |
| 16 | | Geometry correction data record pointer |
| 17 | | Calibration description record pointer |
| 18 | | First calibration data record pointer |
| 19 | | Second calibration data record pointer |
| 20 | | Efficiency pairs record pointer (first record) |
| 21 | | Record number of the first of the three ROI records |
| 22 | | Energy pairs record pointer |
| 23 | | Number of energy pair records |
| 24 | | DLAN1.DLL debug, True = write AN1TRAIL.TXT |
| 25 | | Disable deconvolution of unknown peaks |
| 26 | | True = microcuries, False = becquerels |
| 27 | PERPTR | Laboratory and operator name record pointer |
| 28 | | Maximum record number ever used |
| 29 | | Maximum record number in use |
| 30 | | Number of efficiency pairs records (see Word 20) |
| 31 | SPCTRP | Spectrum record pointer |

| | | |
|-------|--------|-------------------------------------|
| 32 | SPCRCN | Number of records in spectrum |
| 33 | SPCCHN | Number of channels |
| 34 | ABSTCH | Start channel |
| 35–36 | ACQTIM | Acquisition start in DECDAY format |
| 37–40 | ACQT88 | Acquisition start in DECDAY8 format |
| 41 | SEQNUM | Sequence number |
| 42 | MCANUM | Detector ID |
| 43 | SEGNUM | Segment number (always 1) |
| 44 | MCADVT | MCA device type |
| 45 | CHNSRT | Start channel number |
| 46 | RLTMDT | Real time in seconds |
| 48 | LVTMDT | Live time in seconds |
| 50 | | Pointer to MGA records |
| 51 | | Pointer to FRAM records |
| 52 | | Pointer to TRIFID records |
| 53 | | Pointer to NaI records |
| 54 | | Pointer to description record |
| 55–62 | | Reserved |
| 63–64 | | Random summing factor |

Record 2

The second record of the new **.UFO** file contains pointers to other records used in the **.UFO** file and data used by the various analysis programs used to analyze the data.

| <u>Word</u> <u>Number</u> | <u>Local</u> <u>Name</u> | <u>Use</u> |
|------------------------------|-----------------------------|--|
| 1 | I*2 | GEN record pointer |
| 2 | I*2 | CSI 1 record pointer |
| 3 | I*2 | CSI 2 record pointer |
| 4 | I*2 | Library 1 peak records start pointer |
| 5 | I*2 | Number of Library 1 peaks used |
| 6 | I*2 | Number of Library 1 peaks allocated |
| 7 | I*2 | Library 2 peak records start pointer |
| 8 | I*2 | Number of Library 2 peak records used |
| 9 | I*2 | Number of Library 2 peak records allocated |
| 10 | I*2 | Library 3 peak records start pointer |
| 11 | I*2 | Number of Library 3 peak records used |
| 12 | I*2 | Number of Library 3 peak records allocated |
| 13 | I*2 | Unknown peak records start pointer |
| 14 | I*2 | Number of unknown peak records used |

| | | |
|-------|-----|--|
| 15 | I*2 | Number of unknown peak records allocated |
| 16 | I*2 | Library 1 nuclide records start pointer |
| 17 | I*2 | Number of Library 1 nuclide records used |
| 18 | I*2 | Number of Library 1 nuclide records allocated |
| 19 | I*2 | Library 2 nuclide records start pointer |
| 20 | I*2 | Number of Library 2 nuclide records used |
| 21 | I*2 | Number of Library 2 nuclide records allocated |
| 22 | I*2 | Library 3 nuclide records start pointer |
| 23 | I*2 | Number of Library 3 nuclide records used |
| 24 | I*2 | Number of Library 3 nuclide records allocated |
| 25 | | Reserved |
| 28 | I*2 | CSI 3 record pointer |
| 29 | I*2 | CSI 4 record pointer |
| 30 | I*2 | Lowest in-range peak pointer for Lib-1 and this analysis |
| 31 | I*2 | Lowest in-range peak pointer for Lib-2 and this analysis |
| 32 | I*2 | Lowest in-range peak pointer for Lib-3 and this analysis |
| 33 | R*4 | NUPISA |
| 35 | R*4 | EBRDEC |
| 37 | R*4 | IEQVAL |
| 39 | R*4 | IEQDEC |
| 41 | I*2 | NAC shadow for library 1 |
| 43 | I*2 | NAC shadow for library 2 |
| 45 | I*2 | NAC shadow for library 3 |
| 49 | | Pointer to start record 3 |
| 50 | I*2 | Ratio of new nuclide record size to old size. Must be 2 |
| 51 | I*2 | Ratio of new peak record size to old size. Must be 2 |
| 52–58 | | Unused. |
| 59 | R*4 | EBRVAL |
| 61 | R*4 | Total activity |
| 63 | R*4 | Etime |

Other Records

The records not described below are described in the real format spectrum file section.

4.15.2. UFO File CSI Records

Records 1 and 3

UFO File CSI 1 Record

| Word | |
|---------------|-------------------------------------|
| <u>Number</u> | <u>Use</u> |
| 1 | NAA standard table filename |
| 17 | Uranium fission correction filename |
| 33 | Reserved |
| 49 | Reserved |

If long filenames are enabled, then the following 4 sequential records are defined as follows:

| | |
|-----|-------------------------------------|
| 1+2 | NAA standard table filename |
| 3+4 | Uranium fission correction filename |

UFO File CSI 3 Record

| Word | |
|---------------|---------------------------------------|
| <u>Number</u> | <u>Use</u> |
| 1 | MPC Column Title; 64 ASCII characters |
| 33 | Unused |

Records 2 and 4

CSI Records 2 and 4 contain the four filenames that are used in an analysis. Each file specification is stored left-justified in its field as ASCII characters padded with spaces to the end. The contents are:

UFO File CSI 2 Record

| Word | |
|---------------|--------------------|
| <u>Number</u> | <u>Use</u> |
| 1 | RPT filename |
| 32 | MPC Table filename |

If long filenames are enabled, then the following 4 sequential records are defined as follows:

| | |
|-----|--------------------|
| 1+2 | RPT filename |
| 3+4 | MPC Table filename |

UFO File CSI 4 Record (indexed by Word 29 of the second record)

| Word Number | Use |
|----------------|-------------------|
| 1 | UFO filename |
| 32 | Spectrum filename |

If long filenames are enabled, then the following 4 sequential records are defined as follows:

| | |
|-----|-------------------|
| 1+2 | UFO filename |
| 3+4 | Spectrum filename |

4.15.3. UFO File GEN Record

The UFO file GEN record, indexed by Word 1 of the second record, is a special record that contains details about the analysis.

| Word Number | Local Name | Type | Use |
|----------------|---------------|------|-----------------------------------|
| 1–7 | | | Reserved |
| 8 | QA | L*2 | Reserved |
| 9 | IEQ | L*2 | Reserved |
| 10 | EBAR | L*2 | Reserved |
| 11 | RPTLUN | I*2 | Logical unit for report output |
| 12 | UFOLUN | I*2 | Logical unit for UFO file I/O |
| 13 | SPCLUN | I*2 | Logical unit for spectrum I/O |
| 14–15 | | | Reserved |
| 16 | TBLLUN | I*2 | Logical unit for table I/O |
| 17 | LB1LUN | I*2 | Logical unit for Library 1 I/O |
| 18 | LB2LUN | I*2 | Logical unit for Library 2 I/O |
| 19 | LB3LUN | I*2 | Logical unit for Library 3 I/O |
| 20 | LB4LUN | I*2 | Logical unit for Library 4 I/O |
| 21 | PAGENM | I*2 | Current output page number in RPT |
| 22 | LINENM | I*2 | Current output line number in RPT |

| | | | |
|-------|--------|-----|--|
| 23 | NOPIA | I*2 | Number of unknown peaks above sens. cutoff |
| 24 | LIVETM | R*4 | Live time (seconds) value for spectrum |
| 26 | REALTM | R*4 | Real time (seconds) value for spectrum |
| 28 | STTIME | R*4 | Analysis start time |
| 30 | GAMID | I*2 | Type of analysis to perform |
| 31 | GAMVER | C*1 | Analysis version as 8 ASCII characters |
| 35 | GAMNAM | C*1 | ASCII name of analysis program as 8 ASCII characters |
| 39 | ACQTI8 | R*8 | Date and time as double precision DECDAY |
| 43 | DIRFTD | L*2 | True = do directed fit flag |
| 44 | SUM | R*8 | Sum of all channels in the spectrum |
| 49–64 | | | Reserved |

4.15.4. UFO Absorption Description Record

This record is the same as the absorption record in the **.SPC** file (Section 4.11.6).

| <u>Word</u> <u>Number</u> | <u>Local</u> <u>Name</u> | <u>Type</u> | <u>Use</u> |
|------------------------------|-----------------------------|-------------|--|
| 1 | SORPRS | I*2 | Number of pairs (between 1 and 100) |
| 2 | SORINT | I*2 | 1 if geometry record, 2 if internal, 4 if external absorption |
| 3 | SORGEO | I*2 | 1 if cylindrical geometry (internal only), 2 if Marinelli geometry (internal only) |
| 4, 5 | SORDEN | R*4 | Absorber density |
| 6, 7 | SORLNG | R*4 | Absorber length |
| 8–12 | SORDAT | C*10 | Date absorption file created |
| 13–16 | SORTIM | C*10 | Time absorption file created |
| 17–32 | SORNAM | C*32 | Absorber name |
| 33–48 | SREFNAM | C*32 | Reference UFO name |
| 49–64 | SCURNAM | C*32 | Current UFO name |

If long filenames are enabled, then the following 4 sequential records are defined as follows:

| | |
|-----|--------------------|
| 1+2 | Reference UFO name |
| 3+4 | Current UFO name |

4.15.5. UFO Geometry Description Record

Note that this record is the same as the geometry record in the `.SPC` file (Section 4.11.8).

| <u>Word Number</u> | <u>Local Name</u> | <u>Type</u> | <u>Use</u> |
|------------------------|-----------------------|-------------|---|
| 1 | GEOPRS | I*2 | Number of pairs (between 1 and 100) |
| 2 | GEOINT | I*2 | 1 if geometry record, 2 if internal, 4 if external absorption |
| 3 | unused | I*2 | |
| 4–5 | unused | R*4 | |
| 6–7 | unused | R*4 | |
| 8–12 | GEODAT | C*10 | Date geometry file created |
| 13–16 | GEOTIM | C*10 | Time geometry file created |
| 17–32 | unused | C*32 | |
| 33–48 | GREFNAM | C*32 | Reference UFO name |
| 49–64 | GCURNAM | C*32 | Current UFO name |

If long filenames are enabled, then the following 4 sequential records are defined as follows:

| | |
|-----|--------------------|
| 1+2 | Reference UFO name |
| 3+4 | Current UFO name |

4.15.6. UFO File Nuclide Records

4.15.6.1. Old UFO File Nuclide Records

The `.UFO` file nuclide records are 32 words long with two of these records packed in each `.UFO` file 64-word record. These records are allocated in groups of two so that an integer number of 64-word records is always used. The nuclide record contains information about the analysis of a specific nuclide as described below:

| <u>Word Number</u> | <u>Local Name</u> | <u>Type</u> | <u>Use</u> |
|------------------------|-----------------------|-------------|------------------------------------|
| 1 | ISONAM | C*I | Isotope name as 8 ASCII characters |
| 5 | HAFLIF | R*4 | Half-life in days |
| 7 | UNCERT | R*4 | Uncertainty |
| 9 | | I*2 | Reserved |
| 10 | | I*2 | Reserved |
| 11 | PBCVAL | R*4 | PBC correction factor |
| 13 | ISONAA | R*4 | NAA concentration |
| 15 | MPCVAL | R*4 | MPC value |

| 17 | THSREC | I*2 | This record pointer | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
|------------|---|-----|--|------------|--------------------|-------|-----------------------------|-------|------------------------------|-------|----------------------------------|-------|------------------------------|-------|--|-------|--------------------------------------|-------|--|-------|---------------------|-------|--------|-------|---|--------|-------------------------|--------|-----------------------|--------|--------------------------------|--------|----------------|--------|-----------------------------------|--------|--------|
| 18 | ISOFLG | I*2 | Nuclide analysis flag, bits as follows: | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | <table><tr><th><u>Bit</u></th><th><u>Description</u></th></tr><tr><td>Bit 0</td><td>MDA value used for activity</td></tr><tr><td>Bit 1</td><td>Area is below critical level</td></tr><tr><td>Bit 2</td><td>Nuclide must have a report entry</td></tr><tr><td>Bit 3</td><td>No valid peaks were detected</td></tr><tr><td>Bit 4</td><td>Area is between sensitivity and critical level</td></tr><tr><td>Bit 5</td><td>Decay correction time >12 half-lives</td></tr><tr><td>Bit 6</td><td>Decay time during acquisition >12 half-lives</td></tr><tr><td>Bit 7</td><td>Nuclide not present</td></tr><tr><td>Bit 8</td><td>Unused</td></tr><tr><td>Bit 9</td><td>All peaks used in abundance calculation were good</td></tr><tr><td>Bit 10</td><td>Library energy conflict</td></tr><tr><td>Bit 11</td><td>Directed fit activity</td></tr><tr><td>Bit 12</td><td>NAA concentration factor valid</td></tr><tr><td>Bit 13</td><td>NAA calculated</td></tr><tr><td>Bit 14</td><td>Activity calculated using Gamma 3</td></tr><tr><td>Bit 15</td><td>Unused</td></tr></table> | <u>Bit</u> | <u>Description</u> | Bit 0 | MDA value used for activity | Bit 1 | Area is below critical level | Bit 2 | Nuclide must have a report entry | Bit 3 | No valid peaks were detected | Bit 4 | Area is between sensitivity and critical level | Bit 5 | Decay correction time >12 half-lives | Bit 6 | Decay time during acquisition >12 half-lives | Bit 7 | Nuclide not present | Bit 8 | Unused | Bit 9 | All peaks used in abundance calculation were good | Bit 10 | Library energy conflict | Bit 11 | Directed fit activity | Bit 12 | NAA concentration factor valid | Bit 13 | NAA calculated | Bit 14 | Activity calculated using Gamma 3 | Bit 15 | Unused |
| <u>Bit</u> | <u>Description</u> | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Bit 0 | MDA value used for activity | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Bit 1 | Area is below critical level | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Bit 2 | Nuclide must have a report entry | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Bit 3 | No valid peaks were detected | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Bit 4 | Area is between sensitivity and critical level | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Bit 5 | Decay correction time >12 half-lives | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Bit 6 | Decay time during acquisition >12 half-lives | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Bit 7 | Nuclide not present | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Bit 8 | Unused | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Bit 9 | All peaks used in abundance calculation were good | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Bit 10 | Library energy conflict | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Bit 11 | Directed fit activity | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Bit 12 | NAA concentration factor valid | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Bit 13 | NAA calculated | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Bit 14 | Activity calculated using Gamma 3 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Bit 15 | Unused | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 19 | FRPKPT | I*2 | Pointer to first in-range peak for this nuclide | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 20 | | I*2 | Unused | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 21 | | I*2 | Unused | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 22 | | I*2 | Reserved | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 23 | ISOCUR | R*4 | Calculated abundance in becquerels | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 25 | ICRFAC | R*4 | Decay time correction factor | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 27 | ISOERR | R*4 | Isotope counting error | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 29 | TOTERR | R*4 | Total error | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 31 | ISOMDA | R*4 | MDA for this isotope | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |

4.15.6.2. New UFO File Nuclide Records

The first part of the .UFO file nuclide records is 32 words long with two of these records packed in each .UFO file 64-word record (see Section 4.15.1.2). These records are allocated in groups of two so that an integer number of 64-word records is always used. The second part is pointed to by Word 13 below. The nuclide record contains information about the analysis of a specific nuclide as described below:

| <u>Word Number</u> | <u>Local Name</u> | <u>Type</u> | <u>Use</u> |
|------------------------|-----------------------|-------------|---|
| 1 | ISONAM | C*I | Isotope name as 8 ASCII characters |
| 5 | HAFLIF | R*4 | Half-life in days |
| 7 | UNCERT | R*4 | Uncertainty |
| 9 | | I*2 | Reserved |
| 10 | | I*2 | Reserved |
| 11 | PBCVAL | R*4 | PBC correction factor |
| 13 | H2REC | I*2 | Pointer to nuclide record Part 2 |
| 14 | | I*2 | Unused |
| 15 | MPCVAL | R*4 | MPC value |
| 17 | THSREC | I*2 | This record pointer |
| 18 | ISOFLG | I*2 | Nuclide analysis flag, bits as follows: |
| | | <u>Bit</u> | <u>Description</u> |
| | | Bit 0 | MDA value used for activity |
| | | Bit 1 | Area is below critical level |
| | | Bit 2 | Nuclide must have a report entry |
| | | Bit 3 | No valid peaks were detected |
| | | Bit 4 | Area is between sensitivity and critical level |
| | | Bit 5 | Decay correction time >12 half-lives |
| | | Bit 6 | Decay time during acquisition >12 half-lives |
| | | Bit 7 | Nuclide not present |
| | | Bit 8 | Nuclide in table |
| | | Bit 9 | All peaks used in abundance calculation were good |
| | | Bit 10 | Library energy conflict |
| | | Bit 11 | PBC subtracted |
| | | Bit 12 | NAA concentration factor valid |
| | | Bit 13 | Failed key line or fraction test |
| | | Bit 14 | Unused |
| | | Bit 15 | Unused |
| 19 | FRPKPT | I*2 | Pointer to first in-range peak for this nuclide |
| 20 | | I*2 | Unused |
| 21 | | I*2 | Unused |
| 22 | | I*2 | Reserved |
| 23 | ISOCUR | R*4 | Calculated abundance in becquerels |
| 25 | ICRFAC | R*4 | Decay time correction factor |
| 27 | ISOERR | R*4 | Isotope counting error |

| | | | |
|----|--------|-----|----------------------|
| 29 | TOTERR | R*4 | Total error |
| 31 | ISOMDA | R*4 | MDA for this isotope |

Second part:

| <u>Word Number</u> | <u>Local Name</u> | <u>Type</u> | <u>Use</u> |
|------------------------|-----------------------|-------------|--|
| 1 | | | Atomic number and atomic mass of the current isotope |
| 3 | | | Branching ratio to first daughter |
| 5 | | | Branching ratio to second daughter |
| 7 | | | Branching ratio to third daughter |
| 9 | | | Pointer to first parent record |
| 10 | | | Pointer to second parent record |
| 11 | | | Pointer to third parent record |
| 12 | | | Pointer to first daughter record |
| 13 | | | Pointer to second daughter record |
| 14 | | | Pointer to third daughter record |
| 15 | | | Decay during acquisition factor |
| 17 | | | Decay during collection factor |
| 19 | | | Second MDA value (see B30WINS.INI) |
| 21–32 | | | Unused |

4.15.7. UFO File Peak Records

4.15.7.1. Old UFO File Peak Records

The UFO file peak records are 64 words long. These records contain information about the analysis of specific peaks as described below:

| <u>Word Number</u> | <u>Local Name</u> | <u>Type</u> | <u>Use</u> |
|------------------------|-----------------------|-------------|-----------------------------------|
| 1 | PEKENG | R*4 | Library energy of the peak in keV |
| 3 | PGMPRD | R*4 | Gammas per 100 disintegration |
| 5 | PEKCON | R*4 | Peak concentration for NAA |
| 7 | AREA | R*4 | Corrected net area of the peak |
| 9 | BKG | R*4 | Corrected background of the peak |
| 11 | FPAREA | R*4 | First pass net area of the peak |
| 13 | FPBKG | R*4 | First pass background of the peak |

| | | | |
|----|--------|-----|--|
| 15 | CNTR | R*4 | Peak centroid in fractional channels |
| 17 | FWO4M | R*4 | Full-width-.04-maximum in fractional channels |
| 19 | FW10M | R*4 | Full-width-.10-maximum in channels |
| 21 | FWHMOO | R*4 | Full-width-half-maximum in channels |
| 23 | PKUNCT | R*4 | Peak counting uncertainty, as fraction |
| 25 | ISNUM | I*2 | Pointer to nuclide record |
| 26 | ENPRE | I*2 | Pointer to next lower library energy |
| 27 | ENNXT | I*2 | Pointer to next higher library energy |
| 28 | PKARE | I*2 | Pointer to next lower nuclide energy |
| 29 | PKNXT | I*2 | Pointer to next higher nuclide energy |
| 30 | PEK1 | I*2 | Pointer to first in-range peak for this nuclide |
| 31 | FLAG1 | I*2 | Peak flag 1 as described below: |
| | | | <u>Bit</u> <u>Description</u> |
| | | | Bit 0 FW25M/2.2 > 1.2*FWHM |
| | | | Bit 1 FW10M/1.83 > 1.2*FWHM |
| | | | Bit 2 Only 1 entry in library <FWHM |
| | | | Bit 3 2 or more entries in library <FWHM |
| | | | Bit 4 FWHM(actual) < 0.8 FWHM(calc) or > 1.2 FWHM(calc) |
| | | | Bit 5 Peak out of range |
| | | | Bit 6 Net area < 0 |
| | | | Bit 7 Reserved |
| | | | Bit 8 Centroid is too far from library line |
| | | | Bit 9 $HI \leq LO + .8 * FWHM$ |
| | | | Bit 10 Centroid not found |
| | | | Bit 11 Beta peak (obsolete files only) |
| | | | Bit 12 FW25M/2.2 < .8 * FWHM |
| | | | Bit 13 Result of deconvolution |
| | | | Bit 14 Set if error > sensitivity for first in-range peak |
| | | | Bit 15 Unused |
| 32 | FLAG 2 | I*2 | Peak Flag 2 as described below: |
| | | | <u>Bit</u> <u>Description</u> |
| | | | Bit 0 Unused |
| | | | Bit 1 Abundance for this peak higher than average |
| | | | Bit 2 Abundance for this peak lower than average |

| | | | | |
|----|--------|-----|------------|--|
| | | | Bit 3 | Subtraction has occurred |
| | | | Bit 4 | Peak used for MDA calculation |
| | | | Bit 5 | Reserved |
| | | | Bit 6 | Peak used for abundance calculation |
| | | | Bit 7 | Peak used for subtraction |
| | | | Bit 8 | First in-range peak rejected for any reason |
| | | | Bit 9 | Peak itself has error > sensitivity |
| | | | Bit 10 | Duplicate energies in library |
| | | | Bit 11 | Unknown/library multiplet |
| | | | Bit 12 | Peak grossly out of shape |
| | | | Bit 13 | Peak had PBC subtracted |
| | | | Bits 14–15 | Unused |
| 33 | PKCNTS | R*4 | | Number of uncorrected counts in peak |
| 35 | LACTWD | I*2 | | Width of background region below peak |
| 36 | HACTWD | I*2 | | Width of background region above peak |
| 37 | CNTENG | R*4 | | Actual peak centroid in energy |
| 39 | PKLOLM | I*2 | | Peak low limit of integration |
| 40 | PKHILM | I*2 | | Peak high limit of integration |
| 41 | LOAVE | R*4 | | Average background below peak |
| 43 | HIAVE | R*4 | | Average background above peak |
| 45 | PKABUN | R*4 | | Isotopic abundance based on this peak only Or efficiency-corrected area for unknown peaks |
| 47 | PEKMDA | R*4 | | MDA based on this peak only |
| 49 | MPLOLM | I*2 | | Multiplet region low limit of integration |
| 50 | MPHILM | I*2 | | Multiplet region high limit of integration |
| 51 | UFSLOP | R*4 | | Average background slope across peak |
| 53 | PKCONC | R*4 | | |
| 55 | FLAG3 | I*2 | | Peak Flag 3 as described below: |
| | | | <u>Bit</u> | <u>Description</u> |
| | | | Bit 0 | NAC available |
| | | | Bit 1 | Used in concentration calculation |
| | | | Bit 2 | Parabolic background used |
| | | | Bit 3 | Directed fit on this peak |
| | | | Bit 4–15 | Reserved |
| 56 | | I*2 | | Gamma fraction |

| | | | |
|----|--------|-----|---|
| 57 | BKCOEF | R*4 | Parabolic background coefficients (3) |
| 63 | | I*2 | 0 = No other peaks close 1 = Other peaks are too close for reliable deconvolution 2 = The peak area is derived via Gamma3 |

4.15.7.2. New UFO File Peak Records

The UFO file peak records are 64 words long. These records contain information about the analysis of specific peaks as described below:

| <u>Word Number</u> | <u>Local Name</u> | <u>Type</u> | <u>Use</u> |
|--------------------|-------------------|-------------|---|
| 1 | PEKENG | R*4 | Library energy of the peak in keV |
| 3 | PGMPRD | R*4 | Gammas per 100 disintegration |
| 5 | PEKCON | R*4 | Peak concentration for NAA |
| 7 | AREA | R*4 | Corrected net area of the peak |
| 9 | BKG | R*4 | Corrected background of the peak |
| 11 | FPAREA | R*4 | First pass net area of the peak |
| 13 | FPBKG | R*4 | First pass background of the peak |
| 15 | CNTR | R*4 | Peak centroid in fractional channels |
| 17 | FWO4M | R*4 | Full-width-.04-maximum in fractional channels |
| 19 | FW10M | R*4 | Full-width-.10-maximum in channels |
| 21 | FWHMOO | R*4 | Full-width-half-maximum in channels |
| 23 | PKUNCT | R*4 | Peak counting uncertainty, as fraction |
| 25 | ISNUM | I*2 | Pointer to nuclide record |
| 26 | ENPRE | I*2 | Pointer to next lower library energy |
| 27 | ENNXT | I*2 | Pointer to next higher library energy |
| 28 | PKARE | I*2 | Pointer to next lower nuclide energy |
| 29 | PKNXT | I*2 | Pointer to next higher nuclide energy |
| 30 | PEK1 | I*2 | Pointer to first in-range peak for this nuclide |
| 31 | FLAG1 | I*2 | Peak flag 1 as described below: |
| | | <u>Bit</u> | <u>Description</u> |
| | | Bit 0 | FW25M/2.2 > 1.2*FWHM |
| | | Bit 1 | FW10M/1.83 > 1.2*FWHM |
| | | Bit 2 | Only 1 entry in library <FWHM |
| | | Bit 3 | 2 or more entries in library <FWHM |
| | | Bit 4 | FWHM(actual) < 0.8 FWHM(calc) or > 1.2 FWHM(calc) |
| | | Bit 5 | Peak out of range |

| | | | | |
|----|--------|-----|---------------------------------------|--|
| | | | Bit 6 | Net area <0 |
| | | | Bit 7 | Reserved |
| | | | Bit 8 | Centroid is too far from library line |
| | | | Bit 9 | $HI \leq LO + .8 * FWHM$ |
| | | | Bit 10 | Centroid not found |
| | | | Bit 11 | Beta peak (obsolete files only) |
| | | | Bit 12 | $FW25M/2.2 < .8 * FWHM$ |
| | | | Bit 13 | Result of deconvolution |
| | | | Bit 14 | Set if error > sensitivity for first in-range peak |
| | | | Bit 15 | Unused |
| 32 | FLAG 2 | I*2 | Peak Flag 2 as described below: | |
| | | | <u>Bit</u> | <u>Description</u> |
| | | | Bit 0 | Unused |
| | | | Bit 1 | Abundance for this peak higher than average |
| | | | Bit 2 | Abundance for this peak lower than average |
| | | | Bit 3 | Subtraction has occurred |
| | | | Bit 4 | Peak used for MDA calculation |
| | | | Bit 5 | Reserved |
| | | | Bit 6 | Peak used for abundance calculation |
| | | | Bit 7 | Peak used for subtraction |
| | | | Bit 8 | First in-range peak rejected for any reason |
| | | | Bit 9 | Peak itself has error > sensitivity |
| | | | Bit 10 | Duplicate energies in library |
| | | | Bit 11 | Unknown/library multiplet |
| | | | Bit 12 | Peak grossly out of shape |
| | | | Bit 13 | Peak had PBC subtracted |
| | | | Bits 14–15 | Unused |
| 33 | PKCNTS | R*4 | Number of uncorrected counts in peak | |
| 35 | LACTWD | I*2 | Width of background region below peak | |
| 36 | HACTWD | I*2 | Width of background region above peak | |
| 37 | CNTENG | R*4 | Actual peak centroid in energy | |
| 39 | PKLOLM | I*2 | Peak low limit of integration | |
| 40 | PKHILM | I*2 | Peak high limit of integration | |

| | | | |
|--------|--------|-----|---|
| 41 | LOAVE | R*4 | Average background below peak |
| 43 | HIAVE | R*4 | Average background above peak |
| 45 | PKABUN | R*4 | Isotopic abundance based on this peak only or efficiency-corrected area for unknown peaks |
| 47 | PEKMDA | R*4 | MDA based on this peak only |
| 49 | MPLOLM | I*2 | Multiplet region low limit of integration |
| 50 | MPHILM | I*2 | Multiplet region high limit of integration |
| 51 | UFSLOP | R*4 | Average background slope across peak |
| 53 | PKCONC | R*4 | |
| 55 | FLAG3 | I*2 | Peak Flag 3 as described below: |
| | | | <u>Bit</u> <u>Description</u> |
| | | | Bit 0 NAC available |
| | | | Bit 1 Used in concentration calculation |
| | | | Bit 2 Parabolic background used |
| | | | Bit 3 Directed fit on this peak |
| | | | Bit 4–15 Reserved |
| 56 | | I*2 | PKCNUN |
| 57 | BKCOEF | R*4 | Parabolic background coefficients (3) |
| 63 | | I*2 | 0 = No other peaks close 1 = Other peaks are too close for reliable deconvolution 2 = The peak area is derived via Gamma3 |
| 64 | | | Peak flags from library (see lib.h) |
| 65 | MLTPRE | | First peak in multiplet pointer |
| 66 | CURREC | | Current record pointer |
| 67 | LENPRE | | Previous library peak pointer |
| 68 | LENNXT | | Next library peak pointer |
| 69 | G3PRE | | Previous overlapping peak pointer |
| 70 | G3NXT | | Next overlapping peak pointer |
| 71 | FLAGT | | TCC flag |
| | | | <u>Bit</u> <u>Description</u> |
| | | | Bit 1 TCC library |
| 72–73 | MDA2 | | Second MDA value |
| 74–75 | BRUNCT | | Branching ratio uncertainty |
| 76–77 | TCCBR | | Original branching ratio of TCC library peak |
| 78–81 | SUSNAM | | Suspect library nuclide |
| 82–128 | | | Unused |

4.16. Sodium Iodide Unformatted Output (UFO) File Peak Records

4.16.1. UFO Records 1 and 2

The contents of the .UFO file are:

Record 1

| <u>Word</u> <u>Number</u> | <u>Local</u> <u>Name</u> | <u>Use</u> |
|------------------------------|-----------------------------|---|
| 1 | | Must be 1 |
| 2 | | Must be 4096 |
| 3 | | Special application record pointer - 1 |
| 4 | | Special application record pointer - 2 |
| 5 | | Acquisition information record pointer |
| 6 | | Sample description record pointer |
| 7 | | Detector description record pointer |
| 8 | | Special application record pointer |
| 9 | | First analysis parameters record pointer |
| 10 | | Second analysis parameters record pointer |
| 11 | | Third analysis parameters record pointer |
| 12 | | Fourth analysis parameters record pointer |
| 13 | | Absorption correction description record pointer |
| 14 | | Absorption correction data record pointer |
| 15 | | Geometry correction description record pointer |
| 16 | | Geometry correction data record pointer |
| 17 | | Calibration description record pointer |
| 18 | | First calibration data record pointer |
| 19 | | Second calibration data record pointer |
| 20 | | Efficiency pairs record pointer (first record) |
| 21 | | Record number of the first of the two ROI records |
| 22 | | Energy pairs record pointer |
| 23 | | Number of energy pair records |
| 24 | | Reserved |
| 25 | | Disable deconvolution of unknown peaks |
| 26 | | True = microcuries, False = becquerels |
| 27 | PERPTR | Laboratory and operator name record pointer |
| 28 | | Maximum record number ever used |
| 29 | | Maximum record number in use |

| | |
|-------|--|
| 30 | Number of efficiency pairs records (see Word 20) |
| 31–64 | Reserved |

Record 2

The second record of the **.UFO** file contains pointers to other records used in the **.UFO** file and data used by the various analysis programs used to analyze the data.

| Word Number | Local Name | Use |
|----------------|---------------|-----------------------------|
| 1 | | GEN record pointer |
| 2 | | CSI record 1 pointer |
| 3 | | CSI record 2 pointer |
| 4–6 | | Peak record 1 pointer |
| 7–9 | | Peak record 2 pointer |
| 10–12 | | Peak record 3 pointer |
| 13–15 | | unknown peak record pointer |
| 16–18 | | Nuclide record 1 pointer |
| 19–21 | | Nuclide record 2 pointer |
| 22–24 | | Nuclide record 3 pointer |
| 25–27 | | Reserved |
| 28 | | CSI record 3 pointer |
| 29 | | CSI record 4 pointer |
| 30 | LNRGPK | First in-range peak |
| 31 | LNRGP2 | |
| 32 | LNRGP3 | |
| 33–34 | NUPISA | |
| 35–36 | EBRDEC | |
| 37–38 | IEQVAL | |
| 39–40 | IEQDEC | |
| 41–42 | NAC1 | |
| 43–44 | NAC2 | |
| 45–46 | NAC3 | |
| 47–58 | | Reserved |
| 59–60 | EBRVAL | |
| 61–62 | TOTAC | Total activity |
| 63–64 | Etime | Other records |

The remaining record pointers shown in Record 1 above are described in the real format spectrum file section.

4.16.2. UFO File CSI Records

Records 1 and 3

UFO File CSI 1 Record

| Word | |
|---------------|-------------------------------------|
| <u>Number</u> | <u>Use</u> |
| 1 | NAA standard table filename |
| 17 | Uranium fission correction filename |
| 33 | Reserved |
| 49 | Reserved |

If long filenames are enabled, then the following 4 sequential records are defined as follows:

| | |
|-----|-------------------------------------|
| 1+2 | NAA standard table filename |
| 3+4 | Uranium fission correction filename |

UFO File CSI 3 Record

| Word | |
|---------------|---------------------------------------|
| <u>Number</u> | <u>Use</u> |
| 1 | MPC Column Title; 64 ASCII characters |
| 33 | Unused |

Records 2 and 4

CSI Records 2 and 4 contain the four filenames that are used in an analysis. Each file specification is stored left-justified in its field as ASCII characters padded with spaces to the end.

The contents are:

UFO File CSI 2 Record

| Word | |
|---------------|--------------------|
| <u>Number</u> | <u>Use</u> |
| 1 | RPT filename |
| 32 | MPC Table filename |

If long filenames are enabled, then the following 4 sequential records are defined as follows:

| | |
|-----|--------------------|
| 1+2 | RPT filename |
| 3+4 | MPC Table filename |

UFO File CSI 4 Record (indexed by Word 29 of the second record)

| Word | |
|---------------|-------------------|
| <u>Number</u> | <u>Use</u> |
| 1 | UFO filename |
| 32 | Spectrum filename |

If long filenames are enabled, then the following 4 sequential records are defined as follows:

| | |
|-----|-------------------|
| 1+2 | UFO filename |
| 3+4 | Spectrum filename |

4.16.3. UFO File GEN Record

| Word | |
|---------------|-----------------------------------|
| <u>Number</u> | <u>Use</u> |
| 1 | Reserved |
| 2-5 | Reserved |
| 6 | Reserved |
| 7 | Reserved |
| 8 | Reserved |
| 9 | Reserved |
| 10 | Reserved |
| 11 | Logical unit for report output |
| 12 | Logical unit for UFO file I/O |
| 13 | Logical unit for spectrum I/O |
| 14 | Reserved |
| 15 | Reserved |
| 16 | Logical unit for table I/O |
| 17 | Logical unit for library 1 I/O |
| 18 | Logical unit for library 2 I/O |
| 19 | Logical unit For library 3 I/O |
| 20 | Logical unit For library 4 I/O |
| 21 | Current output page number in RPT |
| 22 | Current output line number in RPT |
| 23 | Reserved |

| | |
|-------|--|
| 24 | Live time |
| 26 | Real time |
| 28 | Start time |
| 30 | Type of analysis to perform |
| 31 | Analysis version |
| 35 | Analysis program |
| 39–40 | Acquisition date and time in DECDAY format |
| 41–43 | Reserved |
| 44–45 | Sum of all channels in spectrum |
| 46–47 | Reserved |
| 48–51 | Reserved |
| 52 | Reserved |
| 53 | Reserved |
| 54 | Reserved |
| 55 | Random summing factor |
| 57–58 | Reserved |
| 59–64 | Reserved |

4.16.4. UFO File Nuclide Records

| Word Number | Local Name | Type |
|----------------|--|-----------------------------|
| 1–4 | LIBNAM | |
| 5–6 | HAFLIF | Half-life in days |
| 7–8 | UNCERT | Source uncertainty |
| 9 | NACERR | |
| 10 | THSISO | |
| 11–12 | PBCVAL | |
| 13–14 | ISONAA | |
| 15–16 | MPCVAL | |
| 17 | THSREC | This nuclide record pointer |
| 18 | ISOFLG | Nuclide analysis flag word |
| Bit | Description | |
| 0 | MDA value used for activity | |
| 1 | Area is below critical level | |
| 2 | Nuclide must have a report entry | |
| 3 | No peaks in analysis range | |
| 4 | Area is between sensitivity and critical level | |
| 5 | Decay correction time > 12 half-lives | |
| 6 | Decay during acquisition > 12 half-lives | |
| 7 | Nuclide present | |

| | | | |
|-------|--------|----|---|
| | | 8 | Nuclide in table |
| | | 9 | All peaks used in abundance calculation were good |
| | | 10 | Library energy conflict |
| | | 11 | PBC subtracted |
| | | 12 | NAACON fact available |
| 19 | ISOPEK | | Pointer to first in-range peak for this nuclide |
| 20 | BCKISO | | Pointer to previous isotope nuclide record |
| 21 | FREISO | | Pointer to free isotope nuclide record |
| 22 | NCCNER | | |
| 23–24 | ISOCUR | | Calculated abundance in Bq |
| 25–26 | ICRFAC | | Decay-time correction factor |
| 27–28 | ISOERR | | Isotope counting error |
| 29–30 | TOTERR | | Total error |
| 31–32 | ISOMDA | | MDA for this isotope |

4.16.5. UFO File Peak Records

| Word Number | Local Name | Type |
|----------------|---------------|---|
| 1–2 | PEKENG | Library energy |
| 3–4 | PGMPRD | |
| 5–6 | | Cubic background coefficient GSA change |
| 7–8 | AREA | |
| 9–10 | BKG | |
| 11–12 | AREA1 | (FirstNet) |
| 13–14 | BKG1 | (FirstBackground) |
| 15–16 | CNTR | Centroid in channels |
| 17–18 | FWHM04 | |
| 19–20 | FWHM01 | |
| 21–22 | FWHM00 | |
| 23–24 | PKUNCT | |
| 25 | ISOPTR | Nuclide record (– 1 means this is an unknown peak) |
| 26 | ENPR | Next lower library energy peak record (not necessarily same nuclide) |
| 27 | ENNXT | Next higher library energy peak record (not necessarily same nuclide) |
| 28 | PKPRE | Next lower energy for this nuclide |
| 29 | PKNXT | Next higher energy for this nuclide |
| 30 | PEK1 | Pointer to first in-range peak for this nuclide |

| | | |
|-------|--------|--|
| 31 | FLAG1 | Peak Flag 1 |
| | | <u>Bit</u> <u>Description</u> |
| | | 0 FW25M/2.2 > 1.2 × FWHM |
| | | 1 FW10M/1.83 > 1.2 × FWHM |
| | | 2 Only one entry in library < FWHM |
| | | 3 Two or more entries in library < FWHM |
| | | 4 FWHM(actual) not in range, i.e., !(0.8FWHM(calc) < FWHM(actual) < 1.2FWHM(calc)) |
| | | 5 Peak out of range |
| | | 6 Net area < 0 |
| | | 7 Centroid is too far from library line |
| | | 8 HI ≤ LO + 0.8 × FWHM |
| | | 9 Unused |
| | | 10 Centroid not found |
| | | 11 Beta peak (obsolete) |
| | | 12 FW25M/2.2 < 0.8 × FWHM |
| | | 13 Result of deconvolution |
| | | 14 Error > Sensitivity for first in-range peak |
| 32 | FLAG2 | Peak Flag 2 |
| | | <u>Bit</u> <u>Description</u> |
| | | 0–5 unused |
| | | 6 Used in abundance calculation |
| | | 7 First in-range peak rejected for any reason |
| | | 8 Error > Sensitivity for peak itself |
| | | 9 Peak grossly out of shape |
| | | 10 Unused |
| | | 11 PBC peak |
| 33–34 | PKCNTS | Number of uncorrected counts in peak |
| 35 | | Pointer to next higher present peak |
| 36 | | Number peaks in multiplet region |
| 37–38 | CNTENG | Centroid in energy |
| 39 | PKLOLM | Peak low limit |
| 40 | PKHILM | Peak high limit |
| 41–42 | LOAVE | Average background below peak |
| 43–44 | HIAVE | Average background above peak |
| 45–46 | PKABUN | Isotopic abundance based on this peak only |
| 47–48 | PEKMDA | MDA based on this peak only |
| 49 | MPLOLM | Multiplet low limit |
| 50 | MPHILM | Multiplet high limit |
| 51–52 | UFSLOP | Average background slope above peak (SlopeAbove) |
| 53–54 | PKCONC | (PKCONC) |

| | | |
|-------|--------|--|
| 55 | FLAG3 | Peak Flag 3 |
| | | <u>Bit</u> <u>Description</u> |
| | | 0–3 Unused |
| | | 4 Can use this line to calculate peak contribution |
| | | 5 Subtract peak contribution from this line |
| | | 6 unused |
| | | 7 Use this line to calculate average activity |
| | | 8 Background peak — no activity calculations |
| | | 9 Peak used in fraction limit test, but not for Average Activity |
| | | 10 Peak used in fraction limit test, and, if matched, for Activity |
| | | 11 If nuclide absent, MDA is calculated for this peak energy |
| | | 12 If nuclide absent, LTL is calculated for this peak energy |
| | | 13 If nuclide absent, NET is calculated for this peak energy (and critical level test performed) |
| | | 14 Reserved |
| 56 | | Reserved |
| 57–62 | BKCOEF | Parabolic background coefficients (57,58),(59,60),(61,62) |
| 63 | | Width |
| 64 | | Library peak pointer |

4.17. AlphaVision Output File

AlphaVision stores the analysis results directly in the database and does not create binary output files.

5. ASCII FORMAT FILES

5.1. SPE Spectrum Files

The SPE format files are ASCII text files with several fields before and after the spectrum data. The fields are delimited by fixed keywords beginning with \$ in column 1. The spectrum data is one channel per line. A program can read and use or ignore any fields or keywords not wanted or recognized. Blank lines are ignored. This format is used by several different groups for data interchange.

| <u>Keyword</u> | <u>Content</u> |
|----------------|--|
| \$SPEC_ID: | One line of text describing the data |
| \$SPEC_REM: | Any number of lines containing remarks about the data |
| \$DATE_MEA: | Measurement date in the form mm/dd/yyyy hh:mm:ss |
| \$MEAS_TIM: | Live time and realtime of the spectrum in integer seconds, separated by spaces |
| \$DATA: | The first line contains the channel number of the first channel and the number of channels separated by spaces. The remaining lines contain one channel each of data. |
| \$ROI: | This group contains the regions of interest marked in the spectrum. The first line the number of regions, the following lines contain the start and stop channels for each region. |
| \$ENER_FIT: | This contains the energy calibration factors ($a + b * \text{chn}$) as two real numbers, separated by spaces. |
| \$MCA_CAL: | This contains the number of energy calibration factors on the first line, then the factors on the second line as two numbers, separated by spaces. |
| \$SHAPE_CAL: | This contains the number of FWHM calibration factors on the first line, then the factors on the second line as two numbers, separated by spaces. |

5.2. IEEE Spectrum Format

The IEEE spectrum file consists of a variable number of ASCII lines. Each line is 70 characters long. Each line starts with a 4 character prefix of A004. The prefix is followed by 64 characters of data. All lines end with carriage return and line feed characters. All unused characters in a line are ASCII spaces. The contents of the lines are defined as follows:

| <u>Line</u> | <u>Contents</u> |
|-------------|--|
| 1 | System identification, sub-system identification, analog-to-digital-converter number, segment number, digital offset |
| 2 | Live time, real time, number of channels |
| 3 | Acquisition start time, sample collection time |
| 4 | Energy calibration coefficients |

| | |
|-------|--|
| 5 | Peak full-width-half-maximum (FWHM) calibration coefficients |
| 6 | Sample description — 1 |
| 7 | Sample description — 2 |
| 8 | Sample description — 3 |
| 9 | SAMPLE DESCRIption — 4 |
| 10 | Spare |
| 11–22 | Energy and channel pairs |
| 23–34 | Energy and resolution pairs |
| 35–46 | Energy and efficiency pairs |
| 47–58 | User defined |
| 59 | Spectral data continues to the end |

6. LIBRARY FILES

6.1. Library Files (Germanium)

The germanium library files are created and maintained by the ULI program (A30), the GammaVision program (A66) or the NuclideNavigator programs (B53, C53) and contain the nuclide information needed by the analysis program to analyze germanium spectra. The library file has a header record, nuclide records, and peak records. All types have 64-word records. The extension should be **.LIB**.

6.1.1. Header Record

Record 1

| Word | Local | | | | | | | | | | |
|---------------|----------------------|-------------|---|------------|--------------------|-----------|----------|--------|---------------------|--------|----------------------|
| <u>Number</u> | <u>Name</u> | <u>Type</u> | <u>Use</u> | | | | | | | | |
| 1 | LIBID | I*2 | Must be 4 | | | | | | | | |
| 2 | LIBID2 | I*2 | Each bit represented below: | | | | | | | | |
| | | | <table><tr><th><u>Bit</u></th><th><u>Description</u></th></tr><tr><td>Bits 0–13</td><td>Reserved</td></tr><tr><td>Bit 14</td><td>Sorted library flag</td></tr><tr><td>Bit 15</td><td>Library corrupt flag</td></tr></table> | <u>Bit</u> | <u>Description</u> | Bits 0–13 | Reserved | Bit 14 | Sorted library flag | Bit 15 | Library corrupt flag |
| <u>Bit</u> | <u>Description</u> | | | | | | | | | | |
| Bits 0–13 | Reserved | | | | | | | | | | |
| Bit 14 | Sorted library flag | | | | | | | | | | |
| Bit 15 | Library corrupt flag | | | | | | | | | | |
| 3 | RECLEN | I*2 | Must be 64 | | | | | | | | |
| 4 | RECTOT | I*2 | Number of absolute records in the file | | | | | | | | |
| 5 | | I*2 | Must be 2 | | | | | | | | |
| 6 | TYPREC | I*2 | Pointer to library type name | | | | | | | | |
| 7–8 | | | Unused | | | | | | | | |
| 9 | ISOMAX | I*2 | Number of nuclide records allocated | | | | | | | | |
| 10 | ISOLEN | I*2 | Must be 21 | | | | | | | | |
| 11 | ISOSTR | I*2 | Start of used nuclide records | | | | | | | | |
| 12 | ISOEND | I*2 | End of used nuclide records | | | | | | | | |
| 13 | ISOFRE | I*2 | Start of free nuclide records | | | | | | | | |
| 14 | | I*2 | Must be 1 | | | | | | | | |
| 15 | | I*2 | Must be 1 | | | | | | | | |
| 16 | ISOBAS | I*2 | First 64-word record with nuclide records | | | | | | | | |
| 17 | PKMAX | I*2 | Number of peak records allocated | | | | | | | | |
| 18 | PKLEN | I*2 | Must be 16 | | | | | | | | |
| 19 | PKSTR | I*2 | Start of used peak records | | | | | | | | |
| 20 | PKEND | I*2 | End of used peak records | | | | | | | | |
| 21 | PKFRE | I*2 | Start of free peak records | | | | | | | | |

| | | | |
|-------|--------|-----|---|
| 22 | | I*2 | Must be 2 |
| 23 | | I*2 | Must be 0 |
| 24 | PKBAS | I*2 | First 64-word record with peak records |
| 25–32 | | | Reserved |
| 33 | | C*1 | Creation date and time as 18 ASCII characters. The * character should be ignored if it is not a "1". If it is a "1", it indicates the data is after the year 2000. |
| 43 | | C*1 | Last edit date and time as 18 ASCII characters |
| 52 | NUCUSE | I*2 | Number of used nuclide records |
| 53 | PKUSE | I*2 | Number of used peak records |
| 54–64 | | | Reserved |

6.1.2. Nuclide Records

The nuclide records are each 21 words long. Three of these records are stored on each 64-word record with one word at the end of each 64-word record left unused. Nuclide records are always referred to by a number which is the ordinal number of the record counting from the start of the nuclide records (not a number of a 64-word record).

| <u>Word</u> <u>Number</u> | <u>Local</u> <u>Name</u> | <u>Type</u> | <u>Use</u> |
|------------------------------|-----------------------------|-------------|--|
| 1 | NAME | C*1 | Nuclide name as 8 ASCII characters |
| 5 | HAFLIF | R*4 | Half-life |
| 7 | UNCERT | R*4 | Gammas per disintegration uncertainty |
| 9 | UNI | I*2 | Half-life units (S = 1, M = 2, H = 3, D = 4, or Y = 5) |
| 10 | ISONUM | I*2 | Position in the nuclide chain |
| 11 | | I*2 | Isotope type |
| 12 | | | Nuclide flags |
| | | <u>Bit</u> | <u>Description</u> |
| | | Bit 0 | Thermal neutron activation |
| | | Bit 1 | Fast neutron activation |
| | | Bit 2 | Fission product |
| | | Bit 3 | Naturally occurring isotope |
| | | Bit 4 | Photon reaction |
| | | Bit 5 | Charged particle reaction |
| | | Bit 6 | No MDA calculation |
| | | Bit 7 | Activity not in total |

| | | |
|-------|-----|--|
| 13–16 | | Reserved |
| 17 | I*2 | Pointer to this record |
| 18 | | Reserved |
| 19 | I*2 | Pointer to first peak for this nuclide |
| 20 | I*2 | Back pointer in nuclide chain |
| 21 | I*2 | Fore pointer in nuclide chain |

6.1.3. Peak Records

The peak records are each 16 words long. Four of these records are stored on each 64-word record. Peak records are always referred to by a number which is the ordinal number of the record counting from the start of the peak records (not a number of a 64-word record).

| <u>Word</u> <u>Number</u> | <u>Local</u> <u>Name</u> | <u>Type</u> | <u>Use</u> |
|------------------------------|-----------------------------|-------------|---------------------------------------|
| 1 | KEV | R*4 | Energy in keV |
| 3 | GAMPRD | R*4 | Gammas per disintegration |
| 5 | | I*2 | Peak flags |
| | | | <u>Bit</u> <u>Description</u> |
| | | | Bit 0 Gamma ray |
| | | | Bit 1 X-ray |
| | | | Bit 2 Positron decay |
| | | | Bit 3 Single escape |
| | | | Bit 4 Double escape |
| | | | Bit 5 Key line |
| | | | Bit 6 Not in average |
| 6–8 | | | Unused |
| 9 | PKNUM | I*2 | Position in the sorted peak chain |
| 10 | | I*2 | Pointer to the nuclide record |
| 11–12 | | | Reserved |
| 13 | | I*2 | Back pointer in the sorted peak chain |
| 14 | | I*2 | Fore pointer in the sorted peak chain |
| 15 | | I*2 | Back pointer in the nuclide chain |
| 16 | | I*2 | Fore pointer in the nuclide chain |

6.2. PBC Files (Germanium)

The peak background correction (**.PBC**) files are created and maintained by the GammaVision program (A66) and contain the nuclide information needed by the analysis program to correct for isotopes in the background spectrum. The **.PBC** file has a header record, nuclide records, and peak records. All types have 64-word records. The extension should be **.PBC**.

6.2.1. Header Record

Record 1

| Word | Local | | | | | | | | | | |
|---------------|--------------------|-------------|---|------------|--------------------|-----------|----------|--------|------------------|--------|-------------------|
| <u>Number</u> | <u>Name</u> | <u>Type</u> | <u>Use</u> | | | | | | | | |
| 1 | LIBID | I*2 | Must be 16 | | | | | | | | |
| 2 | LIBID2 | I*2 | Each bit represented below: | | | | | | | | |
| | | | <table><tr><th><u>Bit</u></th><th><u>Description</u></th></tr><tr><td>Bits 0–13</td><td>Reserved</td></tr><tr><td>Bit 14</td><td>Sorted file flag</td></tr><tr><td>Bit 15</td><td>File corrupt flag</td></tr></table> | <u>Bit</u> | <u>Description</u> | Bits 0–13 | Reserved | Bit 14 | Sorted file flag | Bit 15 | File corrupt flag |
| <u>Bit</u> | <u>Description</u> | | | | | | | | | | |
| Bits 0–13 | Reserved | | | | | | | | | | |
| Bit 14 | Sorted file flag | | | | | | | | | | |
| Bit 15 | File corrupt flag | | | | | | | | | | |
| 3 | RECLEN | I*2 | Must be 64 | | | | | | | | |
| 4 | RECTOT | I*2 | Number of absolute records in the file | | | | | | | | |
| 5 | | I*2 | Must be 2 | | | | | | | | |
| 6 | TYPREC | I*2 | Pointer to file type name | | | | | | | | |
| 7–8 | | | Unused | | | | | | | | |
| 9 | ISOMAX | I*2 | Number of nuclide records allocated | | | | | | | | |
| 10 | ISOLEN | I*2 | Must be 21 | | | | | | | | |
| 11 | ISOSTR | I*2 | Start of used nuclide records | | | | | | | | |
| 12 | ISOEND | I*2 | End of used nuclide records | | | | | | | | |
| 13 | ISOFRE | I*2 | Start of free nuclide records | | | | | | | | |
| 14 | | I*2 | Must be 1 | | | | | | | | |
| 15 | | I*2 | Must be 1 | | | | | | | | |
| 16 | ISOBAS | I*2 | First 64-word record with nuclide records | | | | | | | | |
| 17 | PKMAX | I*2 | Number of peak records allocated | | | | | | | | |
| 18 | PKLEN | I*2 | Must be 16 | | | | | | | | |
| 19 | PKSTR | I*2 | Start of used peak records | | | | | | | | |
| 20 | PKEND | I*2 | End of used peak records | | | | | | | | |
| 21 | PKFRE | I*2 | Start of free peak records | | | | | | | | |
| 22 | | I*2 | Must be 2 | | | | | | | | |
| 23 | | I*2 | Must be 0 | | | | | | | | |
| 24 | PKBAS | I*2 | First 64-word record with peak records | | | | | | | | |
| 25–32 | | | Reserved | | | | | | | | |
| 33 | | C*1 | Creation date and time as 18 ASCII characters. The * character should be ignored if it is not a "1". If it is a "1", it indicates the data is after the year 2000. | | | | | | | | |
| 43 | | C*1 | Last edit date and time as 18 ASCII characters | | | | | | | | |

| | | | |
|-------|--------|-----|--------------------------------|
| 52 | NUCUSE | I*2 | Number of used nuclide records |
| 53 | PKUSE | I*2 | Number of used peak records |
| 54–64 | | | Reserved |

6.2.2. Nuclide Records

The nuclide records are each 21 words long. Three of these records are stored on each 64-word record with one word at the end of each 64-word record left unused. Nuclide records are always referred to by a number which is the ordinal number of the record counting from the start of the nuclide records (not a number of a 64-word record).

| Word Number | Local Name | Type | Use |
|----------------|---------------|------|--|
| 1 | NAME | C*1 | Nuclide name as 8 ASCII characters |
| 5–9 | | | Reserved |
| 10 | ISONUM | I*2 | Position in the nuclide chain |
| 11–16 | | | Reserved |
| 17 | | I*2 | Pointer to this record |
| 18 | | | Reserved |
| 19 | | I*2 | Pointer to first peak for this nuclide |
| 20 | | I*2 | Back pointer in nuclide chain |
| 21 | | I*2 | Fore pointer in nuclide chain |

6.2.3. Peak Records

The peak records are each 16 words long. Four of these records are stored on each 64-word record. Peak records are always referred to by a number which is the ordinal number of the record counting from the start of the peak records (not a number of a 64-word record).

| Word Number | Local Name | Type | Use |
|----------------|---------------|------|---------------------------------------|
| 1 | PBCKEV | R*4 | Peak energy in keV |
| 3 | PBCVAL | R*4 | Background peak in counts per second |
| 5–7 | | | Unused |
| 9 | PKNUM | I*2 | Position in the sorted peak chain |
| 10 | | I*2 | Pointer to the nuclide record |
| 11–12 | | | Reserved |
| 13 | | I*2 | Back pointer in the sorted peak chain |
| 14 | | I*2 | Fore pointer in the sorted peak chain |
| 15 | | I*2 | Back pointer in the nuclide chain |
| 16 | | I*2 | Fore pointer in the nuclide chain |

6.3. Library Files in NuclideNavigator

The NuclideNavigator database library files are created and maintained by the NuclideNavigator programs (B53, C53). These files are in Microsoft Access® format. Access can be used to view these files and their structure. NuclideNavigator can also produce GammaVision format libraries as described above.

6.4. Library Files in AlphaVision

AlphaVision (A36-B32) will read NuclideNavigator libraries and convert them to the following format used internally.

| Word Number | Use |
|----------------|---|
| 0 | AlphaVision ID, must be Ox53495641 |
| 1 | Length of each library nuclide field |
| 2 | Number of nuclides (records) in library |
| 3 | Starting record number |
| 4 | Starting peak number |
| 5–14 | Creation date/time |
| 15–24 | Last modification date/time |

6.4.1. Nuclide Records

The nuclide records follow directly after the library header record. The number of records is word 2 of the header.

| Word Number | Use |
|----------------|--|
| 0 | Alpha library entry ID, must be ox53495641 |
| 1–4 | Nuclide name |
| 5–6 | Half-life |
| 7 | Half-life units |
| | 1 seconds |
| | 2 minutes |
| | 3 hours |
| | 4 days |
| | 5 years |
| 8–9 | First peak — energy value |
| 10–11 | First peak — branching ratio |
| 12 | First peak — low-channel ROI value |

| | |
|-------|---|
| 13 | First peak — high-channel ROI value |
| 14 | First peak — starting record number |
| 15 | First peak — starting peak number |
| 16–17 | Second peak — energy value |
| 18–19 | Second peak — branching ratio |
| 20 | Second peak — low-channel ROI value |
| 21 | Second peak — high-channel ROI value |
| 22 | Second peak — starting record number |
| 23 | Second peak — starting peak number |
| 24–25 | Third peak — energy value |
| 26–27 | Third peak — branching ratio |
| 28 | Third peak — low-channel ROI value |
| 29 | Third peak — high-channel ROI value |
| 30 | Third peak — starting record number |
| 31 | Third peak — starting peak number |
| 32–33 | Fourth peak — energy value |
| 34–35 | Fourth peak — branching ratio |
| 36 | Fourth peak — low-channel ROI value |
| 37 | Fourth peak — high-channel ROI value |
| 38 | Fourth peak — starting record number |
| 39 | Fourth peak — starting peak number |
| 40 | Low-channel tracer ROI value |
| 41 | High-channel tracer ROI value |
| 42–45 | Contaminant name |
| 46–47 | Contaminant percent |
| 48–49 | Unused |
| 50–51 | Expected DPM (used in efficiency calibration) |
| 52 | Expected CHN (used in energy calibration) |
| 53–64 | Reserved |

7. DOS SYSTEM FILES

This section describes the files that are used by the system to pass information among programs or to define the hardware system. These are used in DOS, but are not used in the Windows programs.

7.1. Sequence Number File

This file is created by the program SEQNUM (A18) and is always put in directory MCA. It consists of one 64-word direct access record. The values are used to store the spectrum sequence number that is added to the `.SPC` file by CONVERT (A18) or used to generate the default filename in STORE (A18).

The contents are:

| <u>Word Number</u> | <u>Use</u> |
|------------------------|---------------------------------------|
| 1 | Sequence number for MCA #1 Segment 1 |
| 2 | Sequence number for MCA #1 Segment 2 |
| 3 | Sequence number for MCA #1 Segment 3 |
| 4 | Sequence number for MCA #1 Segment 4 |
| . | . |
| . | . |
| . | . |
| . | . |
| 64 | Sequence number for MCA #4 Segment 16 |

7.2. Configuration File

The configuration file is created by the program SETUP (A18) and contains the information about the number of MCAs on the system and the number of segments within those MCAs. The configuration of a system may contain up to 16 MCAs, each having up to 16 segments. The configuration file is always named `\MCA\CONFIG.DAT`. Each record is 256 bytes long.

Record 1

| <u>Word Number</u> | <u>Use</u> |
|------------------------|------------------------|
| 1 | Must be 16 |
| 2 | Number of MCAs defined |
| 3–64 | Unused |

7.2.1. Configuration Records

There is one configuration record for each MCA in the system. The MCA number (as seen by the operator) is one less than the record number of the corresponding configuration record. The contents are:

| <u>Word</u> <u>Number</u> | <u>Type</u> | <u>Use</u> |
|------------------------------|-------------|--|
| 1 | I*2 | Total number of channels in MCA |
| 2 | I*2 | Number of Segments in MCA |
| 3 | C*1 | MCA device name as 4 ASCII characters |
| 5 | I*2 | MCA unit number (0 or 1) |
| 6 | I*2 | MCA type |
| 7 | I*2 | Multiplexer address (0 for none) |
| 8 | | Reserved |
| 9 | | Reserved |
| 10 | | Reserved |
| 12 | I*2 | Start channel number for Segment 1 |
| 13 | I*2 | Number of channels in Segment 1 |
| 14 | | Reserved |
| 15 | | Reserved |
| 17 | I*2 | Start channel number for Segment 2 |
| 18 | I*2 | Number of channels in Segment 2 |
| . | . | . |
| . | . | . |
| . | . | . |
| 84 | | Reserved |
| 85 | I*2 | Start channel number of Segment 16 |
| 87 | I*2 | Number of channels in Segment 16 |
| 88 | I*2 | |
| 89 | | Reserved |
| 90 | C*1 | MCA description as 64 ASCII characters |

7.3. Report File

The report file is created by the program RPT to hold the analysis output report. The extension should be **.RPT**. This file is a formatted, sequential file of variable length records (i.e., it has interrecord delimiters). The DOS commands TYPE or PRINT may be used to read this file.

8. TABLE FILES

8.1. Header Record

These files are mainly used in the DOS programs. Newer formats have been adopted for the Windows programs.

A table file can be created by the programs PBCTABLE or MPCTABLE and contains a table of REAL*4 values associated with a list of nuclides. The table may contain up to four columns of values for each nuclide and it may have up to 300 rows of nuclides and values. A table file has 12-word records with the first five records containing the header information and the rest of the records containing one row each (nuclide name and four values). The following list describes a table file.

Record 1

| Word <u>Number</u> | <u>Use</u> |
|-----------------------|----------------------------|
| 1 | Must be 64 |
| 2 | Reserved |
| 3 | Number of columns in table |
| 4 | Number of rows in table |
| 5–12 | Reserved |

8.2. Title Records

These are Records 2 through 5 and are the titles for Columns 1 through 4, respectively. Title records for unused columns are undefined. They are:

| Word <u>Number</u> | <u>Use</u> |
|-----------------------|-------------------------------------|
| 1 | Column title as 12 ASCII characters |
| 6 | Column units as 12 ASCII characters |

8.3. Nuclide Records

There is one nuclide record for each nuclide defined above. Column values for undefined columns are undefined.

| Word | | |
|---------------|-------------|------------------------------------|
| <u>Number</u> | <u>Type</u> | <u>Use</u> |
| 1 | C*1 | Nuclide name as 8 ASCII characters |
| 5 | R*4 | Column #1 value |
| 7 | R*4 | Column #2 value |
| 9 | R*4 | Column #3 value |
| 11 | R*4 | Column #4 value |

8.3.1. PBC Type

Record 1

| Word | |
|---------------|----------------------------|
| <u>Number</u> | <u>Use</u> |
| 1 | Must be 64 |
| 2 | Reserved |
| 3 | Number of columns in table |
| 4 | Number of rows in table |
| 5–12 | Reserved |

For **.PBC** files the nuclide records are:

| Word | | |
|---------------|-------------|------------------------------------|
| <u>Number</u> | <u>Type</u> | <u>Use</u> |
| 1 | C*1 | Nuclide name as 8 ASCII characters |
| 5 | R*4 | PBC value |
| 7 | R*4 | Unused |
| 9 | R*4 | Unused |
| 11 | R*4 | Unused |